



**So small,  
it can fit in your pocket.**

**So accurate,  
it could go in a lab.**

The new Xprecia Stride™ Coagulation Analyzer  
brings accuracy to the point of care.



### System Features

- ① Test Strip Port
- ② Test Strip Eject Button
- ③ Patient Test
- ④ Recall Results
- ⑤ Settings
- ⑥ Quality Control Test
- ⑦ USB/Micro USB Connector Port
- ⑧ Bar-code Scanner

# System Specifications



<b>System Description</b>	Point-of-care blood coagulation analyzer for PT/INR testing
<b>Calibration</b>	Automatic calibration when a test strip is scanned prior to each test
<b>Sample Type</b>	Fresh capillary (fingerstick) whole blood
<b>Sample Size</b>	Capillary: 6 µL
<b>Interferences</b>	See test strip package insert for details.
<b>Performance Characteristics</b>	Reporting range: 0.8–8.0 INR. Results outside this range will be displayed as out of range.
<b>Memory Recall</b>	640 patient results, 300 liquid quality control results, 300 error messages
<b>Xprecia™ Data Management Software (DMS)</b>	Provides a simple, easy way to transfer data to and from the Xprecia Stride Coagulation Analyzer and a PC. The analyzer connects to a computer via USB to transmit: <ul style="list-style-type: none"> <li>• Operator ID list</li> <li>• Analyzer and operator settings</li> <li>• Analyzer software (firmware) upgrades</li> <li>• Patient and liquid quality control results</li> <li>• Analyzer fault log</li> </ul>
<b>Results Displayed</b>	International normalized ratio (INR), calculated seconds
<b>Electronic Quality Control (EQC)</b>	Two quality-control checks verify the integrity of the test strips. The first control checks the presence of adequate sample and reagents on the test strip, and the second control detects test strip degradation due to exposure to environmental conditions. If either control fails, the system reports an error rather than an incorrect result.
<b>Liquid Quality Control (LQC)</b>	LQC tests help meet regulatory compliance requirements applicable to your facility. Use control solutions to perform quality-control checks on both the analyzer and test strips to ensure they are functioning correctly.

System Dimensions	
<b>Depth</b>	40 mm (1.6 in.)
<b>Height</b>	170 mm (6.7 in.)
<b>Width</b>	70 mm (2.8 in.)
<b>Weight (with batteries)</b>	300 g (10.6 oz)

Xprecia, Xprecia Stride, and all associated marks are trademarks of Siemens Healthcare Diagnostics Inc., or its affiliates. All other trademarks and brands are the property of their respective owners.

Product availability may vary from country to country and is subject to varying regulatory requirements. Please contact your local representative for availability.

## Environmental Requirements

<b>Operating temperature:</b>	15–35°C	59–95°F
20–80% RH, noncondensing at 40°C (104°F)		

<b>Transport and storage:</b>	-20–40°C	-4–104°F
20–85% RH, noncondensing		

**Design life:**  
At least 7,300 cleaning and disinfection cycles

Supported Bar-code Formats	1D	2D
<b>Codabar</b>	•	
<b>Interleaved 2 of 5</b>	•	
<b>Code 39</b>	•	
<b>Code 128</b>	•	
<b>Code 93</b>	•	
<b>Code 49</b>	•	
<b>Data Matrix (ECC200)</b>		•
<b>Aztec</b>		•

## Analyzer Power Source

Disposable alkaline batteries  
Rechargeable nickel metal batteries

**Siemens Healthineers Headquarters**  
Siemens Healthcare GmbH  
Henkestr. 127  
91052 Erlangen  
Germany  
Phone: +49 9131 84 0  
siemens.com/healthineers

**Local Contact Information**  
Siemens Healthcare Diagnostics Inc.  
Point of Care Diagnostics  
2 Edgewater Drive  
Norwood, MA 02062-4637  
USA  
Phone: +1 781-269-3000  
siemens.com/healthineers