


STEP 1

Cleanse wound thoroughly, using local protocols. The skin surrounding the wound should be clean and dry.




STEP 2

Dressing size is dependent on the wound surface area – the dressing pad should extend 1-2cms beyond the wound margin.




STEP 3

Open the sterile packaging and remove the dressing. Separate the paper backing from the AsGUARD® Flex + Non-woven Island Dressing.



STEP 4


Apply the AsGUARD® Flex + Non-woven Island Dressing with central pad covering entire wound bed. Apply gentle pressure around the edge of pad.



STEP 5

Smooth down the edges of the non-woven border ensuring adherence.

DRESSING CHANGE



The clinician-in-charge is responsible for determining the need for dressing changes, dependent on the stage and phase of wound healing and exudate level.

The dressing should be changed when exudate absorption is considered to be extending beyond the wound surface area.

7 days is the maximum period between dressing changes.

To minimise potential damage to fragile skin, AsGUARD® Flex + Non-woven Island Dressing should be eased off using an adhesive remover wipe.

INDICATIONS

AsGUARD® Flex + Non-woven Island Dressings are indicated for use as primary dressings in the management of all types of superficial wounds with light exudate and where water resistance is not required, such as:

- Abrasions and lacerations
- Blisters and minor skin loss
- Skin tears
- Post-operative & post-procedural sites

CONTRAINDICATIONS

- Should not be used if allergies to any of its components, is known.

GENERAL INFORMATION

- If infection is present, a primary wound cleansing dressing e.g. Osmonate® Calcium Alginate dressings or rope, or antimicrobial layer e.g. Zorflex should be used as an alternative.