



Material-friendly rapid cleaning and disinfection wipes impregnated with a low-alcohol solution.

mikrozid® surface wipes

Our Plus

- good material compatibility (e.g. tablets)
- rapid effect
- effective according to EN13727 (bactericidal) in 15 seconds
- effective according to EN16615 (4-field test) in 1 minute
- dermatologically tested
- high quality wipes
- colour and fragrance free
- outstanding cleaning performance

Application areas

Alcohol-based rapid disinfection and cleaning of wipeable surfaces with an increased risk of infection, where short contact times are required and where a material friendly disinfectant is needed, for example:

- Patient and resident nearby surfaces
- Patient and resident treatment areas
- High touch surfaces
- Examination tables and other nearby work surfaces
- Surfaces of medical equipment
- Keyboard or touch pads of sensitive communication equipment, such as smartphones or tablets

Instructions for use

Perform hand hygiene before opening the softpack. Thoroughly wipe the surface with the impregnated cloth and allow to take effect. Make sure to wet surfaces completely and keep them wet for the whole exposure time. Make sure that all visible dirt is removed prior to disinfection. Make sure to close the lid tightly after use.

Not for final disinfection of semicritical and critical medical devices.

If aligned with site policy, this product can be used without gloves.



mikrozid® surface wipes

Product data

Composition:

One softpack contains 120 ready-to-use pre-saturated wipes. 100 g solution contain the following active ingredients: 17,4 g Propan-2-ol, 12,6 g Ethanol (94 % w/w)
Labelling according to Regulation (EC) No. 648/2004: < 5% anionic surfactants.

Chemical / physical data:

Color	colourless
Density	ca. 0,95 g/cm ³
Flash point	26 °C / Method : DIN 51755 Part 1
Form	liquid
pH	3 - 3,6
Viscosity, dynamic	not determined

Special advice

Use disinfectants safely.

Always read the label and product information before use.

The product is highly compatible with metals and synthetics. Please follow carefully any national fire and explosion protection guidelines when using alcohol-based disinfectants. Only treat small surfaces. The wiping range depends on the room temperature and on the structure of the surface to be disinfected. More information available on request.

Information for ordering

Item	Delivery form	Item no.
mikrozid® surface wipes (120 wipes, each wipe 20cm x 20cm)	6/Carton	70002782

Environmental information

schulke manufactures products economically and with advanced, safe and environmentally friendly production processes while at the same time maintaining our high quality standards.