

# MICRO-TOUCH<sup>®</sup> DermaClean<sup>®</sup> Sterile

Performance Series

## Latex, powder-free, sterile, examination glove with excellent grip properties

This sterile glove has excellent grip levels making it ideal for healthcare providers who handle instrumentation, test tubes or infusion equipment in sterile examination (non-surgical) environments. Ansell's natural rubber latex formulation ensures a comfortable, well-fitting glove while the smooth chlorinated inner surface facilitates ease of donning.


### KEY FEATURES AND BENEFITS

- Textured outer surface provides an excellent grip in wet and dry conditions
- Smooth inner surface for easy donning
- Low allergenicity latex reduces risk of contact dermatitis due to absence of thiurams and mercaptobenzothiazoles
- Packaged as a single sterile glove per pouch

Caution: This product contains natural rubber latex which may cause allergic reactions including anaphylactic responses.

### Recommended For

- Sterile procedures outside of the O.R.
- Brief, low-risk examination procedures



## MICRO-TOUCH<sup>®</sup> Everyday, Everywhere Protection

MICRO-TOUCH<sup>®</sup> premium medical examination gloves are scientifically designed to complement every need of medical professionals. Sophisticated design features include Ansell ergonomic technologies that ensure comfort during extended procedures, wet grip innovations that prevent harmful mishandling, and micro-thin films for maximum sensitivity. Caregivers who are fully prepared and comfortably protected by versatile MICRO-TOUCH styles have the peace of mind to focus on what matters most: the health and protection of themselves and their patients.

# MICRO-TOUCH<sup>®</sup> DermaClean<sup>®</sup> Sterile

Performance Series

## PRODUCT INFORMATION

|   |                      |                                    |                                       |
|---|----------------------|------------------------------------|---------------------------------------|
| Material                                | Natural Rubber Latex | Not Made From Natural Rubber Latex | No                                    |
| Color                                   | White                | Cuff Length                        | Standard                              |
| Powder Content                          | Powder-Free          | External Glove Surface             | Fully Textured                        |
| Freedom from Holes (Inspection level I) | 1.5 AQL              | Palm Thickness (mm/mil)            | 0.15 / 5.9                            |
| Finger Thickness (mm/mil)               | 0.17 / 6.7           | Available Sizes                    | S (6.5 - 7), M (7.5 - 8), L (8.5 - 9) |
| Tested For Use With Chemotherapy Drugs  | No                   | Sterile                            | Yes                                   |
| Antistatic                              | Not Tested           | Glove Length (mm/inches)           | 245 / 9.6                             |
| Product Segmentation                    | Performance Series   |                                    |                                       |

## ORDERING INFORMATION

| Size         | S (6.5 - 7) | M (7.5 - 8) | L (8.5 - 9) |
|--------------|-------------|-------------|-------------|
| Product Code | 4480        | 4482        | 4484        |

## PACKAGING AND STORAGE

|                      |  |
|----------------------|--|
| Packaging            | 1 glove /paper pouch 100 gloves per box 4 boxes per carton/case  |
| Shelf Life           | 3 Years  |
| Storage Instructions | Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition. |

## STANDARDS AND CERTIFICATIONS

EN 455 Part 3, EN 420:2003 + A1:2009, ISO 10993-10, ISO 11193, ISO 9001, ISO 13485

Contact your Ansell representative for ordering or more information.

North America  
US Tel: 800 952 9916  
CA Tel: 1-844-494-7854  
Email: insidesalesus@ansell.com

Latin America & Caribbean  
Tel: +52(442) 296 20 50  
Email: cslac@ansell.com  
Brazil Tel: +55-11-3356-3100  
Brazil Email: luvas.medicas@ansell.com

Europe, Middle East & Africa  
Tel: +32 (0) 2 528 74 00  
Email: info@ansell.eu

Asia Pacific  
Tel: +603 8310 6688  
Email: apac.medical@ansell.com

Australia & New Zealand  
Tel: +61 3 9270 7270  
Email: protection@ap.ansell.com