

MICRO-TOUCH[®] DermaClean[®]

Essential Series

Latex, powder-free, examination glove with excellent grip properties

This non-sterile glove has excellent grip levels making it ideal for healthcare providers who handle instrumentation, test tubes or infusion equipment. Ansell's natural rubber latex formulation ensures a comfortable, well-fitting glove while the smooth chlorinated inner surface facilitates ease of donning.


KEY FEATURES AND BENEFITS

- Textured outer surface provides an excellent grip in wet and dry conditions
- Smooth inner surface for easy donning
- Low allergenicity latex reduces risk of contact dermatitis due to absence of thiurams

Caution: This product contains natural rubber latex which may cause allergic reactions including anaphylactic responses.

Recommended For

- Brief, low-risk examination procedures



MICRO-TOUCH[®] Everyday, Everywhere Protection

MICRO-TOUCH[®] premium medical examination gloves are scientifically designed to complement every need of medical professionals. Sophisticated design features include Ansell ergonomic technologies that ensure comfort during extended procedures, wet grip innovations that prevent harmful mishandling, and micro-thin films for maximum sensitivity. Caregivers who are fully prepared and comfortably protected by versatile MICRO-TOUCH styles have the peace of mind to focus on what matters most: the health and protection of themselves and their patients.

MICRO-TOUCH[®] DermaClean[®]

Essential Series

PRODUCT INFORMATION

Material	Natural Rubber Latex	Not Made From Natural Rubber Latex	No
Color	White	Cuff Length	Standard
Powder Content	Powder-Free	External Glove Surface	Fully Textured
Freedom from Holes (Inspection level I)	1.5 AQL	Palm Thickness (mm/mil)	0.14 / 5.5
Finger Thickness (mm/mil)	0.16 / 6.3	Available Sizes	XS (5.5 - 6), S (6.5 - 7), M (7.5 - 8), L (8.5 - 9), XL (9.5 - 10)
Tested For Use With Chemotherapy Drugs	No	Silicone Free	Yes
Sterile	No	Antistatic	Not Tested
Glove Length (mm/inches)	240 / 9.5	Product Segmentation	Essential Series

ORDERING INFORMATION

Size	XS (5.5 - 6)	S (6.5 - 7)	M (7.5 - 8)	L (8.5 - 9)	XL (9.5 - 10)
Product Code	4570	4572	4574	4576	4578

PACKAGING AND STORAGE

Packaging	100 gloves per dispenser 10 dispensers per case 1000 gloves per case *80 gloves per dispenser (Size XL) *800 gloves per case (Size XL)
Shelf Life	3 Years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.

STANDARDS AND CERTIFICATIONS

CA Prop 65, EN 420:2003 + A1:2009, ISO 10993-1, ISO 11193-1:2008, ISO 9001, ISO 10993-5, ISO 13485, ISO 10993-10, ISO 10993-10:2002, Food Type: Rubber

Contact your Ansell representative for ordering or more information.

North America
US Tel: 800 952 9916
CA Tel: 1-844-494-7854
Email: insidesalesus@ansell.com

Latin America & Caribbean
Tel: +52(442) 296 20 50
Email: cslac@ansell.com
Brazil Tel: +55-11-3356-3100
Brazil Email: luvas.medicas@ansell.com

Europe, Middle East & Africa
Tel: +32 (0) 2 528 74 00
Email: info@ansell.eu

Asia Pacific
Tel: +603 8310 6688
Email: apac.medical@ansell.com

Australia & New Zealand
Tel: +61 3 9270 7270
Email: protection@ap.ansell.com