



# Paediatric Cannulas

Salter paediatric cannulas are specifically designed for kids. Kids who like to squirm and aren't very happy with plastic tubing on their cheeks. It has been "kid-tested" for comfort and convenience.



The Smaller Salter-Style Cannula is for oxygen flows up to 3LPM and 6 LPM and one piece anatomical design for secure positioning as well as patient comfort and freedom.



Product Code	Description	QTY
1610	Premature Cannula with 7' (2.13m) supply tubing	Box 50
1611	Neonate Cannula with 7' (2.13m) supply tubing	Box 50
1601	Infant Cannula with 7' (2.13m) supply tubing	Box 50
1615	Intermediate Infant Cannula with 7' (2.13m) supply tubing	Box 50
1602	Paediatric / Small Adult Cannula with 7' (2.13m) supply tubing	Box 50