

# Hot & Cold Gel Packs

## Reusable Hot/Cold Gel Packs

### Usage

- Microwave, Hot Water or Freeze
- To heat, microwave or place in hot water, please follow the instructions on the back of the compress.
- To obtain cold therapy, freeze for at least 1 hour.
- This gel pack can reach a high of 60°C and a low of -12°C
- Maintains therapy for approximately 30 minutes.

### Details

- Size: H125mm x W250mm
- Bag Style: Cloth Feel
- Pack Size: 30 per carton
- Expiry: 3 years from date of manufacture
- No. of uses: No limit
- Note: Please dispose of gel pack if it is seen to be leaking



Order Code	Description
INTGEL3510	REUSABLE HOT/COLD GEL PACK
INTGELDC	CLOTH COVER FOR THE 3510 GEL PACK WHITE