



Conforms to Contours



Self-Adherent



Compress and Supports

CoWrap Cohesive Bandages are constructed from a blend of non-woven fibres and adhesive which adheres to itself, therefore do not require clips. Cohesive bandages conform to body contour movements without bunching or rolling down.

NO TACKY RESIDUE


Does not stick to hair or skin

SELF-ADHERING

Holds in position without the use of a bandage clip, making it easy to re-apply

HAND TEAR

Easy tear bandage, no need for scissors



Suitable for:

Soft tissue support

Securing primary dressings

Compression

Sprains and strains

| Reorder | Size | Colour | Pcs/Unit | Units/Ctn |
|----------|--------------|--------|----------|-----------|
| COB001 | 2.5cm x 4.5m | Tan | 60 | 10 |
| COB051 | 2.5cm x 4.5m | Blue | 60 | 10 |
| COB002 | 5cm x 4.5m | Tan | 36 | 10 |
| COB052 | 5cm x 4.5m | Blue | 36 | 10 |
| COB003 | 7.5cm x 4.5m | Tan | 24 | 10 |
| COB053 | 7.5cm x 4.5m | Blue | 24 | 10 |
| COB004 | 10cm x 4.5m | Tan | 18 | 10 |
| COB054 | 10cm x 4.5m | Blue | 18 | 10 |
| COB005 | 15cm x 4.5m | Tan | 12 | 10 |
| COB055<< | 15cm x 4.5m | Blue | 12 | 10 |

<< Indent Item

