



Comfortable
for Patient



Fast
Absorbency



Gentle
on Skin

Mediam 135 Low Adherent Dressing - Lite consists of a low profile absorbent cotton inner pad encased between two layers of high density polyethylene film.

LOW-ADHERENT

Low adherence to the wound, therefore minimal trauma and pain upon removal

ABSORBENT INNER PAD

Low profile inner pad, wicks excess exudate away from wound bed

CUT TO SIZE


Maintains integrity when cut to size

Suitable for:

Light exuding and small wounds

First aid applications

Minor cuts and grazes



Reorder	Size	Sterility	Pcs/Unit	Units/Ctn
NAL001	7.5cm x 5cm	Sterile	100	8
NAL002	7.5cm x 7.5cm	Sterile	100	8
NAL003	7.5cm x 10cm	Sterile	100	6
NAL004	7.5cm x 20cm	Sterile	100	4

