

Osmocel[®] Si Silicone Hydroporous Foam



Various Size Options



Moisture Vapour Permeable



Thermal Insulation

Osmocel Si Silicone Hydroporous Foam is an all-in-one dressing that protects both skin and wound tissue.

Absorbing excess exudate, this dressing provides the benefits of modern wound healing and encourages less frequent dressing changes.

POLYURETHANE FOAM LAYER

Maintains an optimal moist wound protective layer

ADHESIVE FREE BACKING

Enhanced MVTR that increases possible wear time

EXTENSIVE SIZE RANGE

Catering for extensive wound sizes and body applications

SUPERABSORBENT LAYER*

Superabsorbent layer manages exudate decreasing unnecessary dressing changes



Suitable for:

Superficial wounds with light to moderate exudate

Arterial, diabetic and venous leg ulcers

Skin tears

Pressure injuries, partial thickness burns and surgical wounds



*Bordered range

Osmocel® Si Silicone Hydroporous Foam - Dressing Options



Reorder	Dressing Size	Foam Size	Pcs/Unit	Units/Ctn	Shape
BORDERED					
AWS002	5cm x 12.5cm	2.5cm x 10cm	5	10	
AWS003	7.5cm x 7.5cm	5cm x 5cm	10	10	
AWS004	10cm x 10cm	6.4cm x 6.4cm	10	10	
AWS007	10cm x 20cm	6.5cm x 16.5cm	5	10	
AWS008	10cm x 25cm	6.5cm x 21.5cm	5	10	
AWS009	10cm x 30cm	6.5cm x 26.5cm	5	10	
AWS005	15cm x 15cm	11.5cm x 11.5cm	5	10	
AWS006	15cm x 20cm	11.5cm x 16.5cm	5	10	
LITE BORDERED					
AWS20505	5cm x 5cm	2.5cm x 2.5cm	10	20	
NON-BORDERED					
AWS050	5cm x 5cm		10	10	
AWS051	7.5cm x 7.5cm		10	10	
AWS052	10cm x 10cm		10	10	
AWS053	15cm x 15cm		10	10	
AWS054	15cm x 20cm		10	10	
PETAL BORDERED					
AWS100	7.5cm x 7.5cm	5cm x 5cm	10	10	
AWS101	10cm x 10cm	6.5cm x 6.5cm	10	10	
AWS102	15cm x 15cm	12.5cm x 12.5cm	5	10	
SACRUM BORDERED					
AWS150	18cm x 18cm	13.5cm x 12.5cm	5	10	
AWS151	23cm x 23cm	17.3cm x 16cm	5	10	
SACRUM BUTTERFLY BORDERED					
AWS160	18cm x 18cm	13.8cm x 12.6cm	5	10	
AWS161	23cm x 23cm	19cm x 17.5cm	5	10	
VARI-SITE NON-BORDERED					
AWS180	12cm x 19cm		5	10	
VARI-SITE BORDERED					
AWS182	18.5cm x 24cm	10cm x 13cm	10	10	
AWS183	23cm x 23.2cm	19cm x 20cm	10	10	