



Primary Dressing



Suitable for Fragile Skin



All Exudate Levels

Sofmesh Si Silicone produces an ideal healing environment where the need for atraumatic removal is desired or necessary.

Sofmesh Si Silicone retains moisture at the wound bed, while the perforations allow excess exudate to pass into a secondary dressing, making it suitable for the management of wounds with light to heavy levels of exudate.

## ONE-SIDED ADHERENCE

Easy to apply without adherence to gloves or instruments

## SILICONE CONTACT LAYER

Prevents adherence to the wound, newly formed tissue and surrounding fragile or sensitive skin


## MESH DESIGN

Allows excess exudate to pass into a secondary dressing

## Suitable for:

Treatment of superficial or deep wounds

A protective interface between the wound and adjunct treatments



Reorder	Size	Sterility	Pcs/Unit	Units/Ctn
SOF0575	5cm x 7.5cm	Sterile	10	20
SOF0810	8cm x 10cm	Sterile	10	20
SOF1015	10cm x 15cm	Sterile	10	20
SOF1020	10cm x 20cm	Sterile	10	20
SOF1525	15cm x 25cm	Sterile	10	20

