



Water Proof
Dressing



Conforms
to Body



Highly
Absorbent

SurgiFlex Si Superabsorbent Silicone dressing is a highly absorbent all-in-one island dressing that protects both skin and wound tissue. Absorbing excess exudate, this dressing provides the benefits of modern wound healing and encourages less frequent dressing changes. SurgiFlex Si is ideal for post-surgical procedures where exudate monitoring and management is key.

SILICONE CONTACT LAYER

Suitable for wound types where protection of newly formed tissue is required

SUPERABSORBENT ISLAND

Multi-layer construction to increase exudate absorption

CONFORMABLE AND FLEXIBLE

Spiral embossed pattern provides increased conformability



Suitable for:

Protection of newly formed tissue

Fragile, sensitive or reactive skin

Post-operative wounds and incisions

Reorder	Dressing Size	Padding Size	Sterility	Pcs/Unit	Units/Ctn
SGF6070	6cm x 7cm	3.5cm x 4.5cm	Sterile	10	10
SGF100100	10cm x 10cm	6.5cm x 6.5cm	Sterile	10	10
SGF100200	10cm x 20cm	5cm x 15cm	Sterile	10	10
SGF100250	10cm x 25cm	5cm x 20cm	Sterile	10	10
SGF100300	10cm x 30cm	5cm x 25cm	Sterile	10	10
SGF100350	10cm x 35cm	5cm x 30cm	Sterile	10	10