



Highly Absorbent



Low Linting



X-Ray Detectable

ThreadX Dental Pressure Packs are constructed from cotton with an x-ray detectable thread which is heat bonded to the gauze. The dental pressure packs are highly absorbent with a lengthy tail/tape for ease of removal.

HIGHLY ABSORBENT

Suitable for controlling saliva and bleeding

X-RAY DETECTABLE

X-ray detectable thread for patient safety



Suitable for:

Dental procedures

Prevent and/or control bleeding

Reorder	Size	Pcs	Feature	Thread Colour	Sterility	Pcs/Unit	Units/Ctn
XDP001	Small	1pc	Cotton Tail	Blue	Sterile	100	1
XDP002	Large	1pc	Cotton Tail	Blue	Sterile	100	1
XDP030	Large	2pcs	Cotton Tape	Black	Sterile	100	1