

Hot and Cold Therapy



add-tech pty ltd ABN 57 003 299 572

Unit 2/35 Memorial Ave Ingleburn NSW 2565 PO Box 828 Ingleburn NSW 1890

tel: 02 9618 2577 fax: 02 9618 2588 email: office@addtech.com.au

www.addtech.com.au

addtech 
Innovative Medical Solutions



Hot and Cold therapy

Add-Tech Pty Ltd is an Australian owned medical sales and service company. Add-Tech has been trading since 1988 providing service and importing quality medical products from USA, Europe, Taiwan, Pakistan and China. Our vision is to provide innovative medical solutions to our customers. Add-Tech distributes throughout Australia, New Zealand and the Middle East. Our customers tell us they value our high levels of expertise and dedicated customer service. At Add-Tech we strive to offer quality products, competitive pricing and high levels of customer satisfaction.

The Add-Tech hot and cold therapy product range is manufactured to the highest standard by a leading Chinese manufacturer with over ten year experience manufacturing and exporting throughout the world. All products undergo a stringent quality assurance process both at factory and warehouse level. All products come with a 30 day customer satisfaction guarantee. Our products include, instant hot and cold compresses, reusable hot/cold gel packs, magic warmers hand/body warmers and maternity products.

We thank you for choosing Add-Tech Pty Ltd as a potential trading partner. We look forward to assisting you in the future with all your medical needs.

Product Series:

Instant Cold Compress

Instant Hot Compress

Reusable Hot/Cold Gel Pack

Magic Warmer

Hand/Body Warmer

Wrap

Cold Use:

The cold packs are used to reduce inflammation and swelling by lowering skin temperature. It also helps to relieve pain by numbing injured tissues and reducing muscle spasms. This product is ideal for sprains, strains, contusions, minor burns, toothaches, simple headaches, insect bites and many first aid needs.

Hot Use:

The hot packs are used on injuries after inflammation and swelling has subsided. Hot therapy usually follows cold therapy. These compresses are a good solution to relieve minor muscle and rheumatic pain, joint pain, stiffness and muscle spasms.

Instant Cold Compress



Usage:

With a simple squeeze of the compress break the inner liquid pouch to mix liquid and harmless chemical), the pouch contents become cold immediately and remains cold 15~40 minutes and achieves a low temperature of 25°F.

Characteristic:

Single use cold application without the drippy mess of ice cubes. Needs no refrigeration, convenient and portable, for use whenever and wherever cold therapy is needed. Primary usage is for outdoor activity injuries and first aid.

Item Number	Description	Net Weight (Gross)	Packaging
ICE0506	5"×6"(12.7×15.3cm)	105g(108g)	100pcs/ctn
ICE0609	6"×9"(15.3×25.4cm)	260g(266g)	48pcs/ctn
ICE0506-box	5"×6"(12.7×15.3cm)	105g(119g)	100boxes/ctn
ICE0609-box	6"×9"(15.3×25.4cm)	260g(284g)	48boxes/ctn
ICE0506-3shell	5"×6"(12.7×15.3cm)	315g(355g)	48shells/ctn
ICE0609-box (R)	6"×9"(15.3×25.4cm)	260g(287g)	48boxes/ctn



ICE0609 ICE0506



ICE0609-Box ICE0506-Box



ICE0609-Box(R) ICE0506-3Shell



HOT0609 HOT0506



HOT0609-Box HOT0506-Box



HOT0609-Box(R) HOT0506-3Shell

Instant HOT Compress



Usage:

Simply squeeze the pouch like the instant cold compress. Instant hot compress can last approximately 15~30 minutes and reaches a high temperature of 130°F.

Characteristic:

Product provides heat on it's own after activation. Convenient and portable, for use whenever and wherever hot therapy is needed.

Item Number	Description	Net Weight (Gross)	Packaging
HOT0506	5"×6"(12.7×15.3cm)	95g(98g)	100pcs/ctn
HOT0609	6"×9"(15.3×25.4cm)	218g(224g)	48pcs/ctn
HOT0506-box	5"×6"(12.7×15.3cm)	95g(109g)	100boxes/ctn
HOT0609-box	6"×9"(15.3×25.4cm)	218g(242g)	48boxes/ctn
HOT0506-3shell	5"×6"(12.7×15.3cm)	295g(325g)	48shells/ctn
HOT0609-box (R)	6"×9"(15.3×25.4cm)	218g(245g)	48boxes/ctn

Three simple steps



1) Activate





2) Wrap



3) Apply



Temperature test:



Reusable Hot/Cold Gel Pack



Usage:

To heat, microwave or boil by following instructions on back of the compress. To obtain cold therapy, freeze for at least one hour. It can reach high of 140°F and low of 10°F, and maintains therapy for approximately 30 minutes.

Characteristic:

Can be reused either for cold or hot.
Remains soft and pliable after freezing.
Cleans easily with soap and water.
Different sizes, shapes and colors are available.
Different packaging materials are available, including nylon, PVC, fabric.

Item Number	Description	Net Weight (Gross)	Packaging
GEL0310-Eye	3"×9.6"(7.5×24.5cm)	60g(75g)	100pcs/ctn
GEL0414-Eye	3.7"×14"(9.5×36cm)	120g(150g)	80pcs/ctn
GEL0507	5"×7"(12.7×17.8cm)	140g(114g)	48pcs/ctn
GEL0510-F	5"×10"(12.7×25.4cm)	300g(323g)	30pcs/ctn
GEL0710-P	6.7"×9.8"(17×25cm)	350g(376g)	24pcs/ctn
GEL0812	8"×12"(19×30.8cm)	175g(327g)	24pcs/ctn
GEL0816	8.5"×16"(21.5×40.8cm)	830g(860g)	12pcs/ctn
GEL0507-3Shell	5"×7"(12.7×17.8cm)	420g(463g)	30shells/ctn
GEL0510-Box	5"×10"(12.7×25.4cm)	300g(345g)	30boxes/ctn
GEL0710-P-Box	6.7"×9.8"(17×25cm)	350g(409g)	24boxes/ctn




GEL0507-3Shell



GEL0510-F

Reusable Hot/Cold Gel Pack



GEL0310-Eye



GEL0414-Eye



GEL0507



GEL0710-P



GEL0812



GEL0510-F



GEL0816



Three simple steps


1) Heat Up/Freeze

2) Wrap


3) Apply

Temperature test:

Temperature Curve for GEL0507 (in refrigerator 2 hours)
Room Temperature: 75.2° F



Temperature Curve for GEL0507 (in microwave 20 seconds)
Room Temperature: 75.2° F



Magic Warmer

Magic Warmer



Usage:

Flex the inside metal disc back and forth until crystals form. It can reach temperature of 130°F. To reuse, put the pack in boiling water for several minutes until all crystals are dissolved.

Characteristic:

Beautiful, portable and reusable heat.
Easy to demonstrate.
Different shapes and colors are available.

Item Number	Description	Net Weight (Gross)	Packaging
MAG0304	2.8"×3.5"(7×9cm)	40g(45g)	200pcs/ctn
MAG0305	3.3"×4.5"(8.5×11.5cm)	125g(135g)	60pcs/ctn
MAG0404-Star	3.7"×3.9"(9.5×10cm)	40g(45g)	200pcs/ctn
MAG0405-Heart	4.2"×4.5"(10.6×11.5cm)	60g(70g)	100pcs/ctn
MAG0405	3.5"×5.1"(9×13cm)	160g(173g)	48pcs/ctn
MAG0606	5.5"×5.5"(14×14cm)	190g(208g)	48pcs/ctn
MAG0507	5"×7.5"(13×19cm)	200g(220g)	48pcs/ctn
MAG0715	7"×15"(18.5×38cm)	900g(980g)	48pcs/ctn
MAG0811	7.6"×11"(19.2×28.8cm)	630g(700g)	12pcs/ctn
MAG0816	8.5"×16"(21.5×40.8cm)	780g(840g)	12pcs/ctn
MAG0405-Box	3.5"×5.1"(9×13cm)	160g(183g)	48boxes/ctn
MAG0811-Box	7.6"×11"(19.2×28.8cm)	630g(733g)	12boxes/ctn



MAG0507



MAG0811



MAG0816



MAG0715



MAG0305



MAG0405-Heart



MAG0404-Star



MAG0606





MAG0405



1) Flex metal disc



2) Crystals form



Two simple steps

Hand Warmer/ Body Warmer



WAR0507

Usage:

Open the outer package and shake gently, leaving the internal pack to activate with the exposure to air. It will gradually get warm and once temperature is reached it will last for around 12-15 hours, providing warming for the user. Temperature can reach 124°F ~155°F.

Characteristic:

Product provides long lasting warm to assist user in cold weather. Can be made with adhesive to apply to body or clothing. Ideal for cold weather activities including skiing, camping, sailing, golfing, hiking, ice fishing, snowmobiling or spectator events. Single Use, Disposable.



WAR0204

Item Number	Description	Net Weight (Gross)	Packing
WAR0204	2"×3.9"(5×10cm)	60g(64g)	250pairs/ctn
WAR0507	5"×6.7"(12.5×17cm)	45g(50g)	240pcs/ctn

Wrap with gel pack

Usage:

To heat, microwave or boil the compress by instructions. To obtain cold therapy, freeze the compress for at least one hour. Then place the compress in the bag of the support for usage.

Characteristic:

Used either for cold or hot. Specially designed for parts of body. Remain soft and pliable after freezing. Different sizes and colors are available.



Item Number	Description	Net Weight (Gross)	Packing
WRA0428-Gel	3.7"×11"(9.5×72cm)	535g(575g)	20pcs/ctn
WRA0731-Gel	7.5"×30.7"(19×78cm)	350g(466g)	24pcs/ctn
WRA0738-Gel	7.5"×37.8"(19×96cm)	350g(500g)	24pcs/ctn



WRA0428-Gel



WRA0731-Gel



WRA0738-Gel




Two simple steps 1) Activate



2) Apply



Temperature test:



New



MAG0405-Heel(Infant Heel Warmer)



ICE0513-Perineal(Perineal Cold Pack)