

AUG / SEP / OCT 2024

My Practice Essentials



An EBOS Healthcare Publication
Primary Care Division Phone: 1800 269 534
E&OE. Offers not available to medical wholesalers and
brokers. All offers valid until: 31st October 2024 unless
otherwise specified or while stocks last.

ALL ENCLOSED PRICES ARE GST EXCLUSIVE

ESSENTIAL TO HEALTH

- 10 % GST will apply to all items unless they are marked with GST FREE
- Freight charges apply for orders under \$250 (GST excluded)
- Prices shown are valid from 1st August 2024 until 31st October 2024

Sentry Woundcare @ EBOS:

Osmocel® Si Silicone Hydroporous Foam

Osmocel® Si Silicone Hydroporous Foam is a dressing that protects the skin and wounds. The Osmocel Si promotes wound healing and excess exudate absorption, to decrease the need for frequent dressing changes.

Osmocel Si Silicone Hydroporous Foam is suitable for skin tears, superficial wounds with light exudate, arterial, diabetic, and venous leg ulcers, pressure injuries, partial thickness burns, and surgical wounds.



EBOS Code	Description	UOM	PRICE
10008567	OSMOCEL SI-SILICONE FOAM LITE 5x5cm	Box/10	\$19.08
10015602	OSMOCEL SI-SILICONE BRDRD 7.5x7.5cm	Box/10	\$39.24
10015603	OSMOCEL SI-SILICONE BRDRD 10x10cm	Box/10	\$56.34
10015604	OSMOCEL SI-SILICONE BRDRD 15CMx15cm	Box/5	\$61.66

IMPERVIA® Sterile Absorbent Dressing

IMPERVIA® Sterile Absorbent Dressings are a highly absorbent inner layer with blue impervious backing enclosed by a soft non-woven fabric, making it suitable for primary or secondary dressing, the protection of clothing and linen and moderate to heavy exuding wounds.



EBOS Code	Description	UOM	PRICE
10015464	IMPERVIA® STERILE ABSORBENT DRESSING 10cmX10cm	Box/50	\$31.44
10015465	IMPERVIA® STERILE ABSORBENT DRESSING 10cm X 20cm	Box/25	\$26.22
10015466	IMPERVIA® STERILE ABSORBENT DRESSING 20cm X 20cm	Box/25	\$46.27

AsGUARD® Flex Fabric Fixation Tape, Non-Woven

AsGUARD® Flex Fabric Fixation Tape is a "cut to length" self-adhesive, non-woven tape that is highly comfortable, conforming and breathable.

It is suitable for fixations of dressings and access devices, and acts as a skin protective layer under more aggressive adhesives. Its flexibility also makes this tape ideal and comfortable against joints and difficult contour areas without movement restriction.



EBOS Code	Description	UOM	PRICE
10015523	ASGUARD FLEX NON WOVEN ROLL 2.5cmX10m	Roll	\$3.70
10015524	ASGUARD FLEX NON WOVEN ROLL 5cmX10m	Roll	\$5.28
10015526	ASGUARD FLEX NON WOVEN ROLL 10cmX10m	Roll	\$10.37
10015527	ASGUARD FLEX NON WOVEN ROLL 15cmX10m	Roll	\$14.58

IMPERVIA® Neosorb Superabsorbent Dressing

IMPERVIA® Neosorb Superabsorbent Dressing quickly absorbs and binds large amounts of wound exudate into its 'heart' of Superabsorbent polymers. It is suitable for heavily exuding wounds, exudate management under compression therapy and suitable for protection of clothing and linen.



EBOS Code	Description	UOM	PRICE
10015467	IMPERVIA® NEOSORB SUPERABSORBENT DRESSING 10cm x 10cm	Box/20	\$38.17
10015473	IMPERVIA® NEOSORB SUPERABSORBENT DRESSING 10cm x 20cm	Box/20	\$54.06
10015469	IMPERVIA® NEOSORB SUPERABSORBENT DRESSING 15cm x 20cm	Box/20	\$63.55
10015471	IMPERVIA® NEOSORB SUPERABSORBENT DRESSING 20cm x 30cm	Box/20	\$97.10





Multi-purpose

3M™ Micropore™ S Surgical Tape Pioneering the science of reliable, yet gentle.

Ideal for all patients, including those with sensitive skin.*



Reliable adhesion



Easy to use



Skin performance



Helps reduce risk



Application: **General securement**

- + Blood draws
- + Non-critical tubes
- + Lightweight dressings
- + IV lines and tubing

* Per HRIPT testing.

Australia | Building A, 1 Rivett Road,
North Ryde, NSW 2113, Australia | 1300 363 454

Sponsor | 3M Australia Pty Limited,
1 Rivett Road, North Ryde, NSW 2113

© Solventum 2024. All rights reserved.
Unauthorised use prohibited. Solventum and the
S logo are trademarks of Solventum or its affiliates.

3M™ Micropore™ S Surgical Tape

THE 3M SURGICAL TAPE IS AN ESSENTIAL TOOL FOR HEALTHCARE PROFESSIONALS SEEKING DEPENDABLE ADHESION COMBINED WITH PATIENT COMFORT AND SKIN PROTECTION IN DIVERSE CLINICAL SETTINGS.

Features and Benefits:

- **Reliable Fixation:** Ideal for short-term and long-term (up to 48 hours) applications, ensuring pliable and dependable adherence.
- **Gentle Removal:** Removes cleanly without causing discomfort or disrupting fragile skin layers, supporting patient comfort during tape changes.
- **User-Friendly:** Designed for ease of use, it can be effortlessly repositioned and torn neatly by hand, streamlining application.
- **Skin Integrity:** Helps maintain skin health and integrity, making it suitable for patients with sensitive skin, including infants, children, and the elderly.
- **Non-Irritating:** Formulated to be non-irritating, minimising the risk of adverse skin reactions.
- **Low Hair Pull:** Causes minimal trauma during removal, reducing discomfort associated with hair removal.
- **Versatile Applications:** Suitable for a variety of clinical uses including general securement, blood draws, lightweight dressings, securing IV lines and tubing, and non-critical tubes.



EBOS Code	Description	UOM	PRICE
10021036	3M™ MICROPORE™ S SURGICAL ADHESIVE TAPE 2.5CM X 5M	Box/12	\$89.06
10021037	3M™ MICROPORE™ S SURGICAL TAPE 25MM X 1.37M	Pkt/100	\$212.43
10008722	3M™ MICROPORE™ S SURGICAL TAPE SINGLE PATIENT USE 5CM X 1.3M	Box/50	\$212.46
10008721	3M™ MICROPORE™ S SURGICAL SILICONE-BASED TAPE 5CM X 5M	Box/6	\$86.60



• Freight charges apply for orders under \$250 (GST excluded)
 • Prices shown are valid from 1st August 2024 until 31st October 2024



Good Hand Hygiene



100% AUSTRALIAN OWNED & MADE



Alcohol based hand sanitiser -
 Contains 70% w/w Ethanol - absolute.



Easy to apply and air dries quickly.



Registered therapeutic: AUST R 179333.



Contains two moisturisers, glycerin and dexpanthenol to help hydrate the skin while sanitising.



Alcohol based hand sanitiser helps kill germs fast without water.









Complies with European Standard BS EN 1276 & BS EN 1500.



Available in: 1L, 375ml & 60ml

SANITISING YOUR HANDS CORRECTLY WITH AQIUM GEL



Always read the label.
 Follow the directions for use.

EBOS Code	Description	EACH
10016259	AQIUM ANTI-BACT HAND SANITISER GEL 60ML	\$2.66
10016341	AQIUM ANTI-BACT HAND SANITISER GEL 375ML	\$5.99



Vaccine Fridge with built in UPS Designed and Assembled in Australia! First in-market to have an in-built battery!



Vacc-Safe® Powerguard: Secure Medical Refrigeration with integrated UPS Battery Back-Up

INTRODUCING VACC-SAFE POWERGUARD: AUSTRALIA'S FIRST SECURE MEDICAL REFRIGERATION WITH INTEGRATED UPS BATTERY BACKUP. DESIGNED AND BUILT FOR SAFETY - BECAUSE YOUR CARE DESERVES THE BEST!

SAY GOODBYE TO POWER OUTAGE WORRIES, BIG BULKY UPS BOXES AND HELLO TO SEAMLESS, RELIABLE MEDICAL STORAGE – YOUR VACCINES' NEW BEST FRIEND!



**SIMPLIFIED INSTALLATION
AND SPACE EFFICIENCY.**



**ENSURE VACCINES ARE
KEPT AT THE OPTIMAL
TEMPERATURE AT ALL
TIMES.**



**BATTERY BACK-UP
SUPPLIED AT A
PRE-CALCULATED TIME.**



**COST-EFFECTIVE
SOLUTION**



**SCAN THE QR CODE FOR MORE INFORMATION AND COMPLETE THE
ONLINE FORM TO SUBMIT YOUR EXPRESSION OF INTEREST!**



- 10 % GST will apply to all items unless they are marked with GST FREE
- Freight charges apply for orders under \$250 (GST excluded)
- Prices shown are valid from 1st August 2024 until 31st October 2024

Refrigeration Essentials

Premium Fridge Range

EXCEPTIONAL QUALITY AND TRULY OUTSTANDING FEATURES, AT AN AFFORDABLE PRICE.

- QCPP Compliant
- Independently tested by The Pharmacy Guild of NSW
- Easy to use customised digital controller with High /Low Memory and easy Reset
- Inbuilt USB Data Logger
- Fully adjustable shelving with front and rear upstands
- 24 month warranty
- External access port as standard
- External Alarm contacts supplied as standard



SCAN THE QR CODE FOR MORE INFORMATION AND COMPLETE THE ONLINE FORM TO SUBMIT YOUR EXPRESSION OF INTEREST!

Specifications:

Model	VACC-SAFE PREMIUM 425L	VACC-SAFE PREMIUM 750L	VACC-SAFE PREMIUM 1000L
Dimensions External (W x D x H) mm	690 x 653 x 1920	970 x 808 x 1920	1220 x 808 x 1920
Gross Litres	416	756	1006
Compliance	QCPP and Cold Chain compliant + Pharmacy Guild Tested	QCPP and Cold Chain compliant + Pharmacy Guild Tested	QCPP and Cold Chain compliant + Pharmacy Guild Tested
Warranty	2 Year Warranty	2 Year Warranty	2 Year Warranty
Weight	88	140	220
Operating Range	2 - 8°	2 - 8°	2 - 8°
Shelves	5	10	10
Lock	Safety lock , door self-lock	Safety lock , door self-lock	Safety lock , door self-lock
Alarm	Power failure, Sensor Failure, Door ajar, High and low temp	Power failure, Sensor Failure, Door ajar, High and low temp	Power failure, Sensor Failure, Door ajar, High and low temp
Rated Power (W)	260	320	350
Monitoring	Built-in USB data logger	Built-in USB data logger	Built-in USB data logger
Cooling	Forced air cooling	Forced air cooling	Forced air cooling
Anti-condensation Door	Electrical heating	Electrical heating	Electrical heating



MESI Essentials

MESI mTABLET ECG Diagnostic System

- The most powerful and simultaneous 12-lead acquisition with the University of Glasgow ECG Analysis Algorithm
- Completely wireless and very portable diagnostic device
- It captures and analyses high-quality ECG data in real time that can be shared and accessed for consultation from any device.
- Advanced diagnostic tools and analysis options available on MESI mRECORDS web platform
- Large screen offering portrait and landscape view with various lead layouts

EBOS Code: 11000476



MESI mTABLET ABI Diagnostic System

- PADsense™ algorithm for detection of severe Peripheral Arterial Disease
- Quick and reliable ABI measurement with pulse waveform interpretation

EBOS Code: 11000475



*Category 2 - DIAGNOSTIC PROCEDURES AND INVESTIGATIONS		
11610	Group	D1 - Miscellaneous Diagnostic Procedures And Investigations
	Subgroup	5 - Vascular
Fee: \$69.75 Benefit: 75% = \$52.35 85% = \$59.30		

MESI mTABLET Spirometry Diagnostic System

- Multiple measurement modes and parameter calculation options
- Automatic selection of the best of all repeated maneuvers with the BestBreath™ detection
- Detailed report with the ability to switch between charts and values for a clear interpretation
- Option to compare different measurements on the same screen through MESI mRECORDS
- Pneumotachograph technology with integrated self-calibration for accurate measurements at any time

*Transducers not included with mTABLET Spirometry Diagnostic System

EBOS Code: 11000441



If you are interested in obtaining a quote, please scan the QR Code ▶



- 10 % GST will apply to all items unless they are marked with GST FREE
- Freight charges apply for orders under \$250 (GST excluded)
- Prices shown are valid from 1st August 2024 until 31st October 2024

MESI Essentials

MESI mTablet ECG Module & 4-Port Charging Station

- Glasgow Interpretation Algorithm included
- ECG signal filters (MESI signal enhancement, high-pass, low-pass, mains, myogram)
- Recording speed and sensitivity settings
- Customisable printouts directly from the mTABLET
- Award-winning design and user-friendly interface
- Patient history and trends already on the tablet

EBOS Code: 11000500



MESI mTABLET ABI Diagnostic Module

- PADsense™ algorithm for detection of severe Peripheral Arterial Disease.
- Quick and reliable ABI measurement with pulse waveform interpretation.
- Pulse waveforms and oscilation graphs
- Advanced review and alert

EBOS Code: 11000538



MESI Spirometry Extension Module for mTABLET, Wireless

The MESI Spiro Extension Module allows general practitioners to perform spirometry with an existing MESI mTablet system, offering essential patient monitoring in a wireless manner. The digital Spirometer utilises Pneumotachograph technology with integrated self-calibration for accurate measurements at any time.

EBOS Code: 10010247



If you are interested in obtaining a quote, please scan the QR Code ▶






Ethicon Suture Essentials

MONOCRYL® Plus - Antibacterial

- Physical Material: Poliglecaprone 25 (copolymer of glycolide and epsilon-caprolactone)
- Construction: Monofilament
- Absorbable / Non-absorbable: Absorbable
- Coating (if applicable): N/A
- Colours: Undyed – beige and Dyed – violet
- Available Size Range: 2 through 6/0
- Tensile Strength (Wound Support): Medium Term (3-4 weeks)

MONOCRYL™ Plus Sutures with IRGACARE® MP (triclosan) are effective against the pathogens most commonly associated with surgical site infections.^{1,2}

EBOS Code	Description	UOM	PRICE
10013930	MONOCRYL® Plus Antibacterial Sutures Undyed 70cm 3-0 PS-1	Box/36	\$609.00
10014023	MONOCRYL® Plus Antibacterial Sutures Undyed 45cm 4-0 PS-2	Box/12	\$213.00
10014020	MONOCRYL® Plus Sutures Antibacterial Undyed 45cm 5-0 P-3	Box/12	\$248.00



1. Wang ZX, Jiang CP, Cao Y, Ding DT: Systematic review and meta analysis of triclosan coated sutures for the prevention of surgical site infection. British Journal of Surgery 2013; 154, 89-100
2. Edmiston CE, Daoud FC, Leaper D. Is there an evidence based argument for embracing an antimicrobial (triclosan) coated suture technology to reduce the risk for surgical site infections? A meta analysis. Surgery 2103:154:89-100

VICRYL® Plus - Antibacterial

- Physical Material: Polyglactin 910 (copolymer of glycolide and lactide)
- Construction: Braided (Monofilament in sizes 9/0 and 10/0)
- Absorbable / Non-absorbable: Absorbable
- Coating (if applicable): 50% Polyglactin 370 and 50% Calcium Stearate
- Colours: Undyed beige and Dyed violet
- Available Size Range: 2 through 5/0
- Tensile Strength (Wound Support) Medium Term (4 weeks)

VICRYL™ Plus Sutures with IRGACARE® MP (triclosan) are effective against the pathogens most commonly associated with Surgical site infections.^{1,2}

EBOS Code	Description	UOM	PRICE
10013825	Coated VICRYL* Plus Antibacterial Sutures 70cm 2-0 SH	Box/36	\$279.00
10013833	Coated VICRYL* Plus Antibacterial Sutures 70cm 3-0 FS-1	Box/12	\$335.00
10013838	Coated VICRYL* Plus Antibacterials Sutures 45cm 4-0 P-3	Box/12	\$281.20



1. Wang ZX, Jiang CP, Cao Y, Ding DT: Systematic review and meta analysis of triclosan coated sutures for the prevention of surgical site infection. British Journal of Surgery 2013; 154, 89-100
2. Edmiston CE, Daoud FC, Leaper D. Is there an evidence based argument for embracing an antimicrobial (triclosan) coated suture technology to reduce the risk for surgical site infections? A meta analysis. Surgery 2103:154:89-100

PERMA-HAND® SILK

- Physical Material: Natural protein fibre of raw silk
- Construction: Braided
- Absorbable / Non-absorbable: Non-absorbable
- Coating (if applicable): Bees wax
- Colours: Black and white
- Available Size Range: 5 through 9/0
- Tensile Strength (Wound Support): Medium to Long Term



EBOS Code	Description	UOM	PRICE
10014405	PERMA-HAND* Silk Sutures Black 45cm 2-0 FS 26mm	Box/36	\$83.80
10014404	PERMA-HAND* Silk Sutures Black 45cm 4-0 FS-2	Box/12	\$94.58
10014339	PERMA-HAND* Silk Sutures Black 30cm 6-0 C-3 13mm	Box/12	\$126.31

ETHILON® Nylon Sutures Black

ETHILON Suture is a nonabsorbable, sterile surgical monofilament suture composed of the long-chain aliphatic polymers Nylon 6 and Nylon 6,6. ETHILON Suture is dyed black or green to enhance visibility in tissue.

ETHILON Suture is indicated for use in general soft tissue approximation and/or ligation, including use in cardiovascular, ophthalmic, and neurological procedures.



EBOS Code	Description	UOM	PRICE
10014392	ETHILON* Nylon Sutures Black 45cm 3-0 FS-1 24mm	Box/12	\$89.15
10014391	ETHILON* Nylon Sutures Black 45cm 4-0 FS-2 19mm	Box/12	\$94.86
10014396	ETHILON* Nylon Black 45cm 5-0 FS-2 19mm	Box/12	\$102.84

Many factors can challenge the outcome of your procedures. Bacteria shouldn't be one of them.

From deep tissue layers to the skin, Ethicon products help protect against bacterial infiltration.*

Upgrading to Ethicon's Plus Antibacterial Sutures couldn't be simpler.

BACKGROUND:

- Surgical site infections are the most common and costly healthcare associated infection in Australia occurring at the rate of 36%¹
- The risk factors associated with acquiring a surgical site infection can be multi-factorial, however the risk is increased with the introduction of any foreign body, an example of which are sutures²
- Ethicon Plus antibacterial sutures have been shown in vitro to inhibit bacterial colonisation on the suture for 7 days or more, for protection against the most common organisms associated with surgical site infections (SSIs)^{3,4,5}

BY REDUCING THIS RISK FACTOR:

Ethicon's Plus Sutures have been proven to reduce the risk of SSI by 28%*

- Plus Sutures contains IRGACARE[®] MP (Triclosan) a broad spectrum antibacterial agent .⁶

PDS[™] Plus

Antibacterial
(Polydioxanone)
Suture



MONOCRYL[™] Plus

Antibacterial
(Poliglecaprone 25)
Suture



Coated
VICRYL[™] Plus

Antibacterial
(Polyglactin 910)
Suture



Stratafix[™]


Symmetric PDS[™] Plus
KNOTLESS TISSUE CONTROL DEVICE

Stratafix[™]

Spiral MONOCRYL[™] Plus
KNOTLESS TISSUE CONTROL DEVICE

Stratafix[™]

Spiral PDS[™] Plus
KNOTLESS TISSUE CONTROL DEVICE



References :

1. Russo PL, Stewardson A, Cheng AC, et al Establishing the prevalence of health-care-associated infections in Australian hospitals: protocol for the Comprehensive Healthcare Associated Infection National
2. Edmiston CE, Daoud FC, Leaper D. Is there an evidence-based argument for embracing an antimicrobial (triclosan)-coated suture technology to reduce the risk for surgical-site infections? A meta-analysis. *Surgery* 2013;154:89100.
3. Ming X, Rothenburger S, Yang D. In vitro antibacterial efficacy of Monocryl Plus Antibacterial Suture (poliglecaprone 25 with triclosan). *Surg Infect (Larchmt)*, 2007;8(2):201-207.
4. Rothenburger S, Spangler D, Bhende S, Burkley D. In vitro antimicrobial evaluation of Coated Vicryl Plus Antibacterial Suture (coated polyglactin 910 with triclosan) using zone of inhibition assays. *Surg Infect (Larchmt)*, 2002;3(suppl):S79-S87
5. Ming X, Rothenburger S, Nichols D. In vivo and in vitro antibacterial efficacy of PDS Plus Antibacterial Suture (polydioxanone with triclosan) suture. *Surg Infect (Larchmt)* 2008;9(4):451-458.
6. STRATAFIX[™] Symmetric PDS[™] Plus Knotless Tissue Control Device, Instructions for Use. Somerville, NJ; Ethicon, Inc

Sponsor: Johnson & Johnson Medical Pty Ltd, 1-5 Khartoum Road, North Ryde, NSW EOS : 244017 2303326 26/3/2023

ETHICON
Johnson & Johnson SURGICAL TECHNOLOGIES

When it comes to wound closure have you considered a Topical Skin Adhesive?



A simple to use option that is cost effective and fast.^{2,3,4}

- Addresses cost and convenience concerns by eliminating the need for return visits to remove sutures and any dressings.
- Offers fast closure of small incisions and lacerations.
- Provides patient comfort by providing flexible closure without the pain or anxiety caused by needles.
- Improved cosmesis.⁴

1. Instructions For Use (IFU) Dermabond Mini, Ethicon Inc.
2. Bhende S, Rothenburger S, Spangler DJ, Dito M. In vitro assessment of microbial barrier properties of DERMABOND[™] Topical Skin Adhesive. *Surg Infect* 2002;3(3):251-257.
3. Ethicon information on file 141605-200528
4. Quinn J, Wells G, Sutcliffe T, et al. A randomized trial comparing octyl cyanoacrylate tissue adhesive and sutures in the management of lacerations. *JAMA* 1997;277(19):1527-1530

For complete indications, contraindications, warnings, precautions, & adverse reactions, please reference full package insert.

DERMABOND[™], the Superhero of Topical Skin Adhesives.



Code	Description	Box/12
10014189	DERMABOND [™] MINI	\$249

ETHICON
Johnson & Johnson SURGICAL TECHNOLOGIES



Equipment Essentials:

Examination Couch - 3 Section



Specifications:

Length	
Total	195cm
Head Section	63cm
Middle Section	63cm
Leg Section	63cm
Width	
	71cm
Height Adjustable	46.5cm - 92.5cm
Loading Weight	225kg
Padding thickness	5cm
Foam Density	36kg/m3



EBOS Code	Description	UOM	PRICE
10016248	3 SECTION ELECTRIC COUCH NAVY	Each	\$1,204.04
10016246	3 SECTION ELECTRIC COUCH BLACK	Each	\$1,204.04
10016247	3 SECTION ELECTRIC CHARCOAL GREY	Each	\$1,204.04

2 Section Hi Lo Electric Exam Couch

- Electric Hi Lo Exam Couch operated by foot switch
- Sleek design with rounded frame legs for safety
- Minimum height 50cm to maximum height 90cm in 30 seconds
- Adjustable head section through dual gas operated struts from -20o to +80o
- Head section has face hole and plug
- Adjustable feet for uneven surfaces
- 7.5cm castors for safety and mobility
- Retractable dual locking castors on all 4 corners



EBOS Code	Description	UOM	PRICE
10016244	2 SECTION HI LO ELECTRIC EXAM COUCH NAVY	Each	\$1,204.04
10016243	2 SECTION HI LO ELECTRIC EXAM COUCH GREY	Each	\$1,216.09

Emergency Trolley



- Defibrillator shelf
- Extension table
- 304 Stainless steel IV pole
- Cardiac CPR board
- Waste bin
- External power holder & plug board
- Oxygen tank holder
- 304 Stainless steel wire basket
- Central or seals lock
- Castors with brakes
- Non slip mat
- Luxurious push handle
- Soft closing drawers
- Dimensions: 650 x 480 x 980mm
- Drawer Sizes:
3 x 110mm, 1 x 135mm, 1 x 195mm
- G.W 55kg
- Central locking system through side key



ALL ACCESSORIES
INCLUDED

\$1,372.⁸⁴

EBOS Code: 10011835



Stainless Steel Trolleys

- 50 x 50 x 90 cm



EBOS Code	Description	UOM	PRICE
10011855	STAINLESS STEEL TROLLEY 1 DRAWER	Each	\$426.98
10011857	STAINLESS STEEL TROLLEY 2 DRAWERS	Each	\$564.78
10011859	STAINLESS STEEL TROLLEY 3 DRAWERS	Each	\$676.05



- 10 % GST will apply to all items unless they are marked with GST FREE
- Freight charges apply for orders under \$250 (GST excluded)
- Prices shown are valid from 1st August 2024 until 31st October 2024

Equipment Essentials:

CU Medical Systems® iPAD Automatic External Defibrillator (AED)

The iPAD AED has an output energy of 150 Joules for adults and 50 Joules for paediatric patients with voice prompts for both patient types. The device offers CPR monitoring, patient analysis of heart rhythms and automatic volume adjusting within an LCD display. The display also indicates the status of the device, battery pads, sensing pads expiry date and pads connection status which all increase its ease of use.



GST FREE

\$2,305⁸²

EBOS Code: 10018889

Accessories

EBOS Code	Description	UOM	PRICE
10018872	DEFIBRILLATOR PADS PAEDIATRIC PADS FOR IPAD CU-SP1	Pkt/2	\$195.52
10018887	CU MEDICAL I-PAD SP1 ADULT ELECTRODE PADS	Pkt/2	\$77.13





Equipment Essentials:

ORDISI Examination Lights

State of the art lighting technology and highest light quality, providing the perfect examination conditions. With its long gooseneck arm and swivel bracket feature, easily adjust the light source to the desired position for a better and more accessible examination for diagnosis and treatment.

Ideal for Head, Body, or Gynae Use.

- Swivel fitting - easily moves from left to right
- Extra long gooseneck arm adjust to any position
- High performance LED
- Cold light without infrared,
- LED life 50,000 hours
- Minimal power consumption
- 12 month warranty



EBOS Code	Description	UOM	RRP	PRICE
10012200	LIGHT ORDISI WALL Inc SWIVEL BRACKET	EACH	\$569	\$544.26
10012208	LIGHT ORDISI EXAMINATION MOBILE	EACH	\$569	\$414.86

Welch Allyn® 3.5V Portable Diagnostic Set, PanOptic Basic

The Welch Allyn 3.5V Diagnostic Set includes a PanOptic Basic LED Ophthalmoscope, MacroView Basic LED Otoscope, one Premium Lithium Ion Plus USB Rechargeable power handle and hard case. The Diagnostic set helps keep the tools physicians need fast, organized and at hand to help examine the eyes, nose and throat. The set features long-lasting power handles and LEDs, for more dependable tools to keep physicians running.

The Lithium-Ion Rechargeable power handle is intended for use as a power source for all Welch Allyn standard 3.5V instrument heads.

Hard case for 3.5V diagnostic set organizes and stores the set's components, creating greater ease for clinicians.



\$1,453.⁶³

EBOS Code: 10009021

WelchAllyn®

- 10 % GST will apply to all items unless they are marked with GST FREE
- Freight charges apply for orders under \$250 (GST excluded)
- Prices shown are valid from 1st August 2024 until 31st October 2024

Diagnostic Essentials

Xprecia Stride™ Coagulation Analyser

Primary care, urgent care, and hospitals demand fast, reliable PT/INR testing for oral anticoagulation therapy (OAT). That's why EBOS Healthcare is introducing the Xprecia Stride™ Coagulation Analyser, a truly handheld device that tests PT/INR with lab accuracy.

The Xprecia Stride analyser is ergonomically designed, with a large, intuitive touchscreen to enhance the user experience, simplify work flow, and improve overall testing accuracy. Together, these features make the testing process seamless and help minimise errors.

\$350

EBOS Code: 10008181



**Only
\$3.20 a
strip!**



**SAVE
on your
strips!**

EBOS Code	Size	UOM	PRICE
10008180	XPRECIA STRIDE PT/INR TEST STRIPS (UBI)	Box/100	\$282.75



Sage®

Introducing our **NEW** Sterile Instrument Dispensers




NEW



EBOS Code	Description	UOM	PRICE
10022665	SAGE® BLUNT DRESSING FORCEPS STERILE SINGLE USE 12.5CM	Box/50	\$71.15
10022667	SAGE® ADSON DRESSING FORCEPS NO TEETH STERILE 12CM	Box/60	\$152.67
10022666	SAGE® ADSON DRESSING FORCEPS 1:2 TEETH STERILE 12CM	Box/60	\$148.99
10022670	SAGE® MOSQUITO HALSTED STRAIGHT FORCEPS STERILE 12.5CM	Box/25	\$87.05
10022669	SAGE® MCINDOE FORCEPS STERILE 15CM	Box/30	\$132.08
10022671	SAGE® MOSQUITO HALSTED CURVED FORCEPS STERILE 12.5CM	Box/20	\$70.97
10022674	SAGE® MAYO HEGAR NEEDLE HOLDER 15CM	Box/20	\$78.35
10022675	SAGE® IRIS PRECISION SCISSORS SHARP BLADE STERILE 11.5CM	Box/50	\$129.74
10022676	SAGE® SHARP BLUNT TIP SCISSORS STERILE 12.5CM	Box/35	\$80.86

Quality, Precision and Convenience in Every Pack





Flaminal®

For a variety of wound types.

Due to it's unique composition

Flaminal® for acute and chronic wounds or wounds at risk of infection.

Exclusively available with EBOS Healthcare

Flaminal® is a **versatile, dependable and unique** dressing designed to simplify wound care. It provides clinical benefits in just one powerful product:



Suitable for many wounds and wound conditions:

NECROTIC WOUNDS

Necrotic material can delay healing and promote infection.



INFECTED WOUNDS

Bacteria may impede healing, exacerbate pain and cause malodour.

Under medical supervision only (ref 14)



SLOUGHY WOUNDS

Yellow slough and necrotic material can delay healing and promote infection. Risk of maceration in exuding wounds.



GRANULATING WOUNDS

Granulation tissue can be fragile and wound edges need protecting.



Flaminal^{hydro}

Indicated for **low to moderately** exuding wounds



Flaminal^{forte}

Indicated for **moderate to highly** exuding wounds

EBOS Code	Description	PRICE Each
10019908	Flaminal Forte 50g	\$22.90
10019909	Flaminal Hydro 50g	\$22.90



Flamigel®

Flamigel® is a safe and effective treatment for superficial wounds.

- Sunburn
- Superficial burn
- Cuts
- Grazes
- Post Laser Therapy
- Radiodermatitis



\$14.90

EBOS Code: 10019903


*Reference: Thomas, S. Wound management and dressings. The pharmaceutical press London. 1990, Chapter 2: Functions of a wound dressing, p9-19.

4-Step wound care with Flamigel®

Flamigel®'s hydro-active colloid formulation helps absorb excess moisture for wet wounds or rehydrates dry wounds. As a result, Flamigel® ensures optimal moisture balance that promotes **40% faster healing***.



1 Cleaning Wash your hands and clean the wound thoroughly. Rinse the wound under running water or use a wound cleanser.




Treatment

- ✓ Apply a layer of Flamigel® 1-2 times per day.
- ✓ For sunburn, apply Flamigel® 2-6 times per day.



2 Drying Gently dry the wound with a clean dry towel.



Cover

If necessary, cover the wound with a plaster or dressing and change daily.

If the wound becomes swollen, reddened or painful, or doesn't start to heal within a few days, seek medical advice.



Surface Cleaning and Disinfecting Essentials

TUFFIE 5 Cleaning and Disinfecting Wipes

The TUFFIE 5 Universal Sanitising, Cleaning and Disinfecting Wipes are powerful alcohol and chlorine-free wipes that are ideal for patient care environment where cleanliness and eradication of bacteria is required. The unique formulation, combined with the mechanical action when using the wipe, aids the removal of organic matter and improves infection prevention and control.



EBOS Code	Description	Tub/225
901SW225	TUFFIE 5 UNIVERSAL WIPES	\$18.95

TUFFIE Detergent Wipes

The TUFFIE Detergent wipes are designed with an improving cleaning performance with non-streak formula, faster drying time and active odour control, making the Tuffie wipes ideal for use on mattresses, beds, commodes, door handles, hoists, trolleys and other general equipment



EBOS Code	Description	Tub/225
901DW225	TUFFIE DETERGENT WIPES	\$17.95

Vernacare

Patient Wipes Essentials

CONTI®

POST TOILETING

Introducing...

Conti® Post Toileting Wet Wipes

A more sustainable cleaning option for assisted intimate hygiene.

Flushable and Maceratable.



EBOS Code	Description	Each
10007666	Conti Cleansing Wet Wipes, Pack of 12 x 50 - Pack of 50. Wipe size 22 x 17cm, Fragrance Free	\$48.80
10021703	Conti Cleansing Wet Wipes, Pack of 27 x 24 - Pack of 24. Wipe size 22 x 17cm, Fragrance Free	\$110.00

- 10 % GST will apply to all items unless they are marked with GST FREE
- Freight charges apply for orders under \$250 (GST excluded)
- Prices shown are valid from 1st August 2024 until 31st October 2024

Vaccine Essentials

AREXVY
(RESPIRATORY SYNCYTIAL VIRUS VACCINE, ADJUVANTED)

10021865	AREXVY RESPIRATORY SYNCYTIAL VIRUS VACCINE, ADJUVANTED	EACH	\$309
----------	--	------	--------------

Vaxchora
(Cholera Vaccine, Live, Oral)

10021836	VAXCHORA CHOLERA VACCINE (ORAL) SACHET	EACH	\$179
----------	--	------	--------------

SHINGRIX
(HERPES ZOSTER VACCINE RECOMBINANT, ADJUVANTED)

10010916	SHINGRIX VIAL + AS01B VIAL AU	EACH	\$268.51
----------	-------------------------------	------	-----------------

havrix Adult
Inactivated hepatitis A virus vaccine

10010911	HAVRIX ADULT PFS 1 ML X 1	EACH	\$51.33
----------	---------------------------	------	----------------

twinrix Adult
combined hepatitis A and hepatitis B vaccine

10010905	TWINRIX ADULT PFS 1 ML X 1	EACH	\$60.35
----------	----------------------------	------	----------------

Trumenba
Meningococcal Group B Vaccine

10011995	TRUMENBA MENINGOCOCCAL B VACCINE	EACH	\$123.69
----------	----------------------------------	------	-----------------

havrix Child
Inactivated hepatitis A virus vaccine

10010912	HAVRIX JNR PFS 1 ML X 1	EACH	\$33.96
----------	-------------------------	------	----------------

twinrix Child
combined hepatitis A and hepatitis B vaccine

10010904	TWINRIX JUNIOR PFS	EACH	\$41.11
----------	--------------------	------	----------------

Adacel
Pertussis Vaccine Acellular Combined with Diphtheria and Tetanus Toxoids (adsorbed)

10019822	ADACEL 0.5ML PFS	EACH	\$28.03
----------	------------------	------	----------------

BEXSERO
Meningococcal Group B Vaccine

10010986	BEXSERO PFS 0.5 ML X 1	EACH	\$121
----------	------------------------	------	--------------

STAMARIL
Yellow Fever Vaccine (Live), Stabilised

10019824	STAMARIL YELLOW FEVER VACCINE	EACH	\$99.85
----------	-------------------------------	------	----------------

Prevnar 20
Pneumococcal 20-valent Conjugate Vaccine

10008262	PREVENAR 20 VAC 1X0.5ML PFS	EACH	\$146.70
----------	-----------------------------	------	-----------------

IPOL
Inactivated poliovirus vaccine (IPV)

10017248	IPOL (SALK-VERO) 0.5 ML	EACH	\$42.68
----------	-------------------------	------	----------------

MSD M-M-R
(Measles, Mumps, & Rubella Virus Vaccine Live)

10017243	MMR-II VIAL X1 WITH PFS	EACH	\$3.64
----------	-------------------------	------	---------------

boostrix
Combined diphtheria-tetanus-acellular pertussis (dTpa) vaccine

10010987	BOOSTRIX PFS 0.5 ML X1	EACH	\$59.35
----------	------------------------	------	----------------

TYPHIM Vi
Salmonella Typhi Vi Polysaccharide Vaccine

10019823	TYPHIM VI 0.5 ML	EACH	\$59.72
----------	------------------	------	----------------

AVAXIM
Inactivated Hepatitis A Vaccine

10019825	AVAXIM HEP A 0.5ML AGE 2+	EACH	\$44.25
----------	---------------------------	------	----------------

engerix b
hepatitis B surface antigen recombinant (yeast)

10010900	ENGERIX-B ADULT PFS 1 ML X1	EACH	\$28.10
----------	-----------------------------	------	----------------

Nimenrix
Meningococcal polysaccharide groups A, C, W-135 and Y conjugate vaccine

10010903	NIMENRIX ACWY DOSE PFS DILUENT	EACH	\$60.42
----------	--------------------------------	------	----------------

Priorix
Measles, mumps, rubella vaccine, live

10010983	PRIO.VIAL +DIL.PRTC AU	EACH	\$22.31
----------	------------------------	------	----------------

H-B-VAX II
Hepatitis B vaccine (recombinant)

10017239	HB VAX II ADULT 10MCG/1ML X1 PFS	EACH	\$20.94
----------	----------------------------------	------	----------------

Vivotif
Typhoid Vaccine Live Oral Ty21a

10008740	VIVOTIF ORAL CAPSULES TYPHOID VACCINE	EACH	\$54.15
----------	---------------------------------------	------	----------------

MENVEO
Meningococcal (Groups A, C, W-135 and Y) Oligosaccharide CRM197 Conjugate Vaccine

10010915	MENVEO VIAL VIAL 0.5 ML X1	EACH	\$60.80
----------	----------------------------	------	----------------



HARTMANN
+

Zetuvit® Plus Silicone Border

For moderately to highly exuding wounds.

- +** Effective
- +** Patient-friendly
- +** Versatile

Designed
to be your
first choice



Hartmann Zetuvit® Plus Silicone Dressing

Zetuvit® Plus Silicone is a soft silicone superabsorbent dressing designed to effectively manage moderate to high levels of exudate in acute and chronic wounds. By locking away bacteria and exudate, it supports better healing conditions and helps prevent complications. Its outer layer is water repellent, and can be worn for up to 7 days (depending on the condition of the wound and the clinician's instructions). Its silicone contact layer is gentle on the skin and allows for atraumatic removal.



EBOS Code	Description	Pkt/10
10008198	HARTMANN ZETUVIT® PLUS SILICONE DRESSING, 8CM X 8CM	\$78.38

Zetuvit Plus Silicone Border Dressing

Zetuvit® Plus Silicone Border is a soft silicone superabsorbent dressing designed to effectively manage moderate to high levels of exudate in acute and chronic wounds. By locking away bacteria and exudate, it supports better healing conditions and helps prevent complications. With its showerproof and secure adhesion, it stays in place for up to 7 days (depending on the condition of the wound and the clinician's instructions). Its silicone contact layer and border is gentle on the skin and allows for atraumatic removal.



EBOS Code	Description	Pkt/10
10015333	ZETUVIT PLUS SILICONE BORDER DRESSING 10CM X 10CM, WHITE	\$58.24
10015336	ZETUVIT PLUS SILICONE BORDER DRESSING 15X25CM	\$233.13





Your Primary Care Team

Western Australia

- **Western Regional Sales Manager - Glen Charles**
gcharles@ebosgroup.com.au
- **Key Account Manager - Deena Doyle**
ddoyle@ebosgroup.com.au
- **Key Account Manager - Joanna Humphries**
jhumphries@ebosgroup.com.au

Queensland

- **Key Account Manager - Patricio Hernandez**
phernandez@ebosgroup.com.au
- **Key Account Manager - Shaun Callander**
scallander@ebosgroup.com.au
- **Key Account Manager - Leanne Carey**
leanne.carey@ebosgroup.com.au

South Australia

- **Key Account Manager - Malcolm Wainwright**
malcolm.wainwright@ebosgroup.com
- **Key Account Manager - Soula Pica**
Soula.Pica@ebosgroup.com

Victoria

- **National Corporate Accounts Manager - Michael Mazengarb**
mmazengarb@ebosgroup.com.au
- **customerservice@ebosgroup.com.au**

New South Wales

- **Key Account Manager - Carol Wellard**
cwellard@ebosgroup.com.au
- **Key Account Manager- Sarah Tiriaki**
stiriaki@ebosgroup.com.au
- **Key Account Manager - Helen Koveos**
hkoveos@ebosgroup.com.au

Place an order by:

- Calling us on **1800 269 534**
- Sending us a Fax on **1800 810 257**
- Sending us an Email
customerservice@ebosgroup.com.au
- Visiting us online
www.eboshealthcare.com.au

EBOS Group Australia Pty Ltd
ABN 38 125 401 247

Availability of products and pricing subject to change without notice.

Free freight on all orders over \$250 (GST Excluded) – excluding equipment indicated with this icon , which will incur additional freight charges. Please call customer service for a quotation.

PRICES in catalogue exclude 10% GST. 10% will apply to all items listed unless marked GST FREE. Valid until stocks last and finishes on 31/10/2024.

Supporting our community for
100 Years