

# Basin Medical Mona Lisa<sup>®</sup> Copper IUDs

- BARIUM SULPHATE FOR RADIOLOGICAL PROOFING.
- MR CONDITIONAL.



## Mona Lisa CU375 and Mona Lisa CU375 SL

This model is identical to the long-term, tried and trusted best seller "Multiload CU 375/375 SL".

Insertion requires a one-hand technique only with preloading.

This Mona Lisa is available in 2 versions:

### Mona Lisa<sup>®</sup> Cu375

- Copper surface: 375 mm<sup>2</sup>
- In-situ-time: 5 years
- One-handed insertion technique

### Mona Lisa<sup>®</sup> Cu375 SL

- Identical characteristics to the Mona Lisa Cu375
- Shorter stem length (29.4 mm)
- Ideally suited for a shorter uterus length

## Mona Lisa CUT 380A QL

This model is an improved version of the "Gyne-T".

The classic T-shaped IUD is the most commonly used IUD model in the world, recommended even by the WHO. Both IUDs feature a secure fit in the uterus, the lowest Pearl-Index and a minimal expulsion rate.

### Mona Lisa<sup>®</sup> CuT 380A QL

- Copper surface: 380 mm<sup>2</sup>
- In-situ-time: 10 years
- IUD with the lowest Pearl Index
- Easy and quick loading



EBOS Code	Supplier Code	Product	Width	Length	Maximum Insertion Duration	Recommended Uterus Length
10022753	BMCU375	MONA LISA <sup>®</sup> CU375	19.5mm	35.2mm	5 Years	6-9cm
10022754	BMCU375-SL	MONA LISA <sup>®</sup> CU375 SL	19.5mm	29.4mm	5 Years	5-8cm
10022752	BMCUT380A-QL	MONA LISA <sup>®</sup> CUT 380A QL	31.5mm	35.9mm	10Years	6-9cm



FOR MORE INFORMATION CONTACT YOUR ACCOUNT MANAGER  
ON **FREEPHONE 1800 269 534** OR VISIT [www.ebosonline.com.au](http://www.ebosonline.com.au)



ESSENTIAL TO HEALTH