

INTRODUCING...

Mona Lisa[®]

Cu375 / 375SL Copper IUDs

FROM **BASIN**  **MEDICAL**

Cu375 / 375 SL

This model is identical to the long-term tried and trusted bestseller "Multiload Cu 375/375 SL".

Insertion requires a one-hand technique only without preloading.

Available in two different versions:

Cu375

With a copper surface of 375 mm² this version guarantees 5 years of safe contraception.

Unique features:

- ✓ *Mona Lisa[®] Cu375: 19.5 mm wide and 35.2 mm long*
- ✓ *The only Australian MR Conditional Copper IUD*
- ✓ *One-hand technique for insertion*
- ✓ *Low expulsion rate due to the flexible IUD-body*

Cu375 SL

This version is identical to the Mona Lisa Cu 375, only differing in the dimension.

Unique features:

- ✓ *Mona Lisa[®] Cu375 SL: 19.5 mm wide and 29.4 mm long. Stem length is 5mm shorter than the Mona Lisa Cu375*
- ✓ *Suited for patients with shorter uterus length*



BASIN  **MEDICAL**

For further information on these products visit the Basin Medical website www.basinmedical.com.au or email sales@basinmedical.com.au