

MIC* & MIC-KEY*
Enterostomy Tubes



INNOVATIVE PRODUCTS G, J AND GJ-TUBES WITH ENFit™ CONNECTION

DIGESTIVE HEALTH PRODUCT CATALOGUE



INNOVATIVE PRODUCTS - G, J AND GJ-TUBES WITH ENFit™ CONNECTION

DIGESTIVE HEALTH PRODUCT CATALOGUE

TABLE OF CONTENTS

INNOVATIVE PRODUCTS FOR G, J AND GJ-TUBES - INTRODUCTION 4

GASTRIC ACCESS FEEDING TUBES

| | |
|--|---|
| MIC-KEY* G-Tube | 6 |
| MIC* G-Tube | 7 |
| MIC* BOLUS G-Tube | 7 |
| MIC* PEG Kit | 8 |
| HALYARD* Introducer Kit for G-Tube | 9 |

JEJUNAL ACCESS FEEDING TUBES

| | |
|-----------------------------|----|
| MIC-KEY* J-Tube | 9 |
| MIC* J-Tube | 10 |
| MIC* Jejunostomy Tube | 10 |

GASTRO-JEJUNAL ACCESS FEEDING TUBES

| | |
|--|----|
| MIC-KEY* GJ-Tube | 11 |
| MIC* GJ-Tube (Endoscopic/Fluoroscopic Placement) | 12 |
| MIC* GJ-Tube (Surgical Placement) | 12 |
| MIC* Gastro-Enteric Tube | 13 |
| HALYARD* Introducer Kit for J-/GJ-Tube | 13 |

ACCESSORIES

| | |
|---|----|
| MIC-KEY* Extension Sets | 14 |
| MIC* Jejunostomy Feeding Tube Accessories | 14 |
| MIC* Extension Sets | 15 |
| MIC* PEG Feeding Connector Sets | 15 |

INNOVATIVE PRODUCTS FOR G, J AND GJ-TUBES - INTRODUCTION



Halyard Health, a leader in the enteral feeding market, offers a wide variety of innovative, high-quality enteral feeding tubes, placement kits and accessories uniquely designed for delivering nutrition and medication to paediatric and adult patients.

As the pioneer in the development of the first balloon retained tube designed for gastrostomy feeding, Halyard continues to manufacture tubes and accessories with the highest standards of workmanship at our state-of-the-art manufacturing facilities.

Our products provide consistent performance and are designed for safety and convenience. For this reason our MIC* and MIC-KEY* feeding tubes are now equipped with the ENFit™ connector. This connector is designed to connect only to other ENFit™ compatible devices, so that only enteral fluids are given down the feeding tube.

We stand behind our products with knowledgeable sales and customer service professionals. Today, our MIC*/MIC-KEY* brand of Enteral Feeding Tubes is well recognised and accepted among healthcare professionals for intermediate and long-term feeding applications.

MIC*/MIC-KEY*

Enterostomy Tubes

INNOVATIVE PRODUCT DESIGNS

MIC*/MIC-KEY* Enteral Feeding Tubes and Accessories include a wide variety of innovative product designs for gastric, jejunal and gastro-jejunal access: from the lock-and-key design of the MIC-KEY* Feeding Tube extension set connection mechanism, to the one-piece MIC-KEY* GJ Feeding Tube that allows simultaneous gastric decompression and jejunal feeding.

MOST MIC*/MIC-KEY* ENTEROSTOMY TUBES INCLUDE:

- High-clarity, medical-grade silicone construction, designed for visibility and drapability;
- Ventilated SECUR-LOK* ring, the external retention mechanism designed to allow air to circulate around the stoma site and reduce tube pressure 360° from the sides of the stoma wall;
- An inflatable internal retention balloon, which reduces patient trauma experienced with obturated devices;
- ENFit™ feeding port and medication port compliant to the ISO80369-3 standard.

THE RIGHT CHOICE

Helping your patients face the use of a gastrostomy tube, possibly for a lifetime, can be a challenge. Halyard's variety of quality feeding tubes and accessories helps you make the right choice for your patient. By selecting MIC*/MIC-KEY* Enterostomy Tubes, you may be assured that you have chosen innovative products that consistently perform and are well recognised and accepted among healthcare professionals.



GASTRIC ACCESS FEEDING TUBES

HALYARD* Gastric Access Feeding Tubes are designed for the delivery of enteral nutrition to the stomach and are available in low-profile and conventional, standard-length designs.

HALYARD* MIC-KEY* G-TUBE

The **MIC-KEY* G-TUBE** is a step up from conventional gastrostomy feeding tubes. The most trusted gastrostomy feeding tube, the MIC-KEY* G Feeding Tube is unobtrusive and easy to conceal, making it an ideal feeding tube for virtually all individuals.

- Medical-Grade Silicone Construction
- Low-Profile Design
- Tapered Distal Tip Matching Tube Size
- Inflatable Silicone Internal Retention Balloon
- Not formulated with natural rubber latex nor DEHP
- ETO sterilised
- With online tutorial for patients and caregivers
- Recessed a-traumatic tip
- High performing Anti-Reflux Valve
- SECUR-LOK* Extension Set Connector Mechanism with ENFit™ Feeding Port
- Radiopaque Striping
- 1 Each/Packaging




HALYARD* MIC-KEY* G-TUBES

| Size (Outer Dia.) (X) | Stoma Length (Y) | Product Code |
|-----------------------|------------------|--------------|
| 12 Fr | From 0.8 to 4.0 | 8140-X-Y |
| 14 Fr | From 0.8 to 5.0 | 8140-X-Y |
| 16 Fr | From 0.8 to 5.0 | 8140-X-Y |
| 18 Fr | From 0.8 to 5.0 | 8140-X-Y |
| 20 Fr | From 0.8 to 5.0 | 8140-X-Y |
| 24 Fr | From 1.5 to 5.0 | 8140-X-Y |

Available stoma length (cm): 0.8 /1.0 /1.2 /1.5 /1.7 /2.0 /2.3 /2.5 /2.7 /3.0 /3.5 /4.0 /4.5 /5.0

Kit includes:

- 1- MIC-KEY* Low-Profile G-tube
- 1- MIC-KEY* ENFit™ Extension Set with Right Angle Connector and double male ENFit™ connector, 30cm length
- 1- MIC-KEY* Bolus Extension Set with Straight Connector and male ENFit™ connector, 30cm length
- 1- 6ml Luer-Slip™ Syringe
- 1- 35ml female ENFit™ syringe
- 4- Gauze Pads




HALYARD* MIC* G-TUBE

The first of its kind, the **MIC* G TUBE** is indicated for patients requiring continuous enteral feeding.

- Medical-Grade Silicone Construction
- Recessed a-traumatic Tip (except Fr 12)
- Inflatable Silicone Internal Retention Balloon (7-10ml, except LV-versions: 3-5ml)
- SECUR-LOK* External Retention Ring
- Not formulated with natural rubber latex nor DEHP
- Gamma sterilised
- Dual Ports: Male ENFit™ Feeding and Medication Port
- Tapered Distal Tip Matching Tube Size
- Radiopaque Striping
- 1 Each/Packaging

HALYARD* MIC* G-TUBES

| Size (Outer Dia.) | Balloon Volume | Product Code |
|-------------------|----------------|--------------|
| 12 Fr | 3-5ml Balloon | 8100-12LV |
| 14 Fr | 3-5ml Balloon | 8100-14LV |
| 16 Fr | 3-5ml Balloon | 8100-16LV |
| 18 Fr | 7-10ml Balloon | 8100-18 |
| 20 Fr | 7-10ml Balloon | 8100-20 |
| 22 Fr | 7-10ml Balloon | 8100-22 |
| 24 Fr | 7-10ml Balloon | 8100-24 |
| 26 Fr | 7-10ml Balloon | 8100-26 |
| 28 Fr | 7-10ml Balloon | 8100-28 |
| 30 Fr | 7-10ml Balloon | 8100-30 |



HALYARD* MIC* BOLUS G-TUBE

The **MIC* BOLUS G FEEDING TUBE** is utilised for patients who are on a bolus feeding regimen.

- Medical-Grade Silicone Construction
- Recessed a-traumatic Tip (except 12 Fr)
- Inflatable Silicone Internal Retention Balloon
- SECUR-LOK* External Retention Ring
- Male ENFit™ Feeding Port
- Tapered Distal Tip matching tube size
- Radiopaque Striping
- Not formulated with natural rubber latex nor DEHP
- Gamma sterilised
- 1 Each/Packaging

HALYARD* MIC* BOLUS G-TUBES

| Size (Outer Dia.) | Balloon Volume | Product Code |
|-------------------|----------------|--------------|
| 12 Fr | 3-5ml Balloon | 8110-12LV |
| 14 Fr | 3-5ml Balloon | 8110-14LV |
| 16 Fr | 3-5ml Balloon | 8110-16LV |
| 18 Fr | 7-10ml Balloon | 8110-18 |
| 20 Fr | 7-10ml Balloon | 8110-20 |
| 22 Fr | 7-10ml Balloon | 8110-22 |
| 24 Fr | 7-10ml Balloon | 8110-24 |



HALYARD* MIC* PERCUTANEOUS ENDOSCOPIC GASTROSTOMY (PEG) KIT

The **MIC* PEG KIT** is designed for placement convenience. Available in PULL and PUSH version.

- Medical-Grade Silicone Construction
- SECUR External Retention Ring
- Male ENFit™ Feeding and Medication Port
- Not formulated with Natural Rubber Latex nor DEHP
- Radiopaque Stripe and full opaque Bumper
- Tubing Clamp
- Traction removable (except Fr 14)
- ETO Sterilised
- 1 Each/Packaging

HALYARD* MIC* PEG KITS

| Size (Outer Dia.) | Description | Code |
|-------------------|---------------------------|---------|
| 14 Fr | PUSH Method | 8630-14 |
| 20 Fr | PUSH Method | 8630-20 |
| 24 Fr | PUSH Method | 8630-24 |
| 14 Fr | PULL Method | 9640-14 |
| 20 Fr | PULL Method | 9640-20 |
| 24 Fr | PULL Method | 9640-24 |
| 14 Fr | PULL Method (limited kit) | 8644-14 |
| 20 Fr | PULL Method (limited kit) | 8644-20 |
| 24 Fr | PULL Method (limited Kit) | 8644-24 |
| 14 Fr | PULL Method | 8645-14 |
| 20 Fr | PULL Method | 8645-20 |
| 24 Fr | PULL Method | 8645-24 |

Series with item#8644-xx contains only products indicated with †

PUSH PEG Kit Includes:

- 1- Traction Removable PUSH PEG Tube
- 1- 5 Fr Retrieval Snare (1.8mm X 260cm)
- 1- ENFit™ Connector
- 1- SECUR-LOK* Retention Ring
- 1- 19ga X 1.5in. Filter Needle
- 1- 25ga X 1.5in. Needle
- 1- Fenestrated Drape with Adhesive Tabs
- 1- 12ml Syringe
- 1- Scalpel with #11 Blade
- 1- Stainless Steel Introducer Needle
- 1- Guidewire
- 1- Tubing Clamp
- 1- Water-Soluble Lubricant
- 4- 4in. X 4in. Gauze Sponges
- 2- 2in. X 2in. Slit Gauze Sponges
- 1- Suture Strand

PULL PEG Kit Includes:

- 1- Traction Removable PULL PEG Tube†
- 1- 5 Fr Retrieval Snare (1.8mm X 260cm) (not available in item#8645-xx)
- 1- ENFit™ Connector†
- 1- SECUR-LOK* Retention Ring†
- 1- 19ga X 1.5in. Filter Needle
- 1- 25ga X 1.5in. Needle
- 1- Fenestrated Drape with Adhesive Tabs
- 1- 12ml Syringe
- 1- Scalpel with #11 Blade†
- 1- Introducer Needle†
- 1- Looped Placement Wire†
- 1- Tubing Clamp†
- 1- Water-Soluble Lubricant
- 4- 4in. X 4in. Gauze Sponges
- 2- 2in. x 2 in. Slit Gauze Sponges
- 1- Suture Strand



Kit includes:

- 1- Gastrointestinal Anchor Set with 4 Pre-loaded SAFE-T-PEXY* T-Fasteners
- 1- Haemostat
- 1- Safety Introducer Needle
- 1- J-tipped Guide wire
- 1- Safety Scalpel #11
- 1- Serial Dilator with Peel-away Sheath
- 2- Syringes 12ml
- 1- Stoma Measuring Device

Feeding tube not included in kit.

HALYARD* INTRODUCER KIT FOR G-TUBE PLACEMENT

Kits include tools and devices to facilitate primary endoscopic, fluoroscopic or laparoscopic placement of balloon retained gastrostomy tubes such as MIC* G-tubes and MIC-KEY* G-tubes.

- ETO Sterilised
- 1 Each/Packaging

HALYARD* INTRODUCER KITS FOR G-TUBE PLACEMENT

| Description | Product Code |
|--|--------------|
| Introducer Kit for G-Tube placement, with 16 Fr Dilator (for use with MIC* G-Tube Fr 12 and MIC-KEY* G-Tube Fr 12) | 98430 |
| Introducer Kit for G-Tube placement, with 18 Fr Dilator (for use with MIC* G-Tube Fr 14 and MIC-KEY* G-Tube Fr 14) | 98431 |
| Introducer Kit for G-Tube placement, with 20 Fr Dilator (for use with MIC* G-Tube Fr 16 and MIC-KEY* G-Tube Fr 16) | 98432 |
| Introducer Kit for G-Tube placement, with 22 Fr Dilator (for use with MIC* G-Tube Fr 18 and MIC-KEY* G-Tube Fr 18) | 98433 |
| Introducer Kit for G-Tube placement, with 24 Fr Dilator (for use with MIC* G-Tube Fr 20 and MIC-KEY* G-Tube Fr 20) | 98434 |

JEJUNAL ACCESS FEEDING TUBES

HALYARD* Jejunal Access Feeding Tubes are designed for the delivery of enteral nutrition to the small bowel and are available in low-profile and conventional designs.



Kit Includes:

- 1- MIC-KEY* Low-Profile J-tube
- 1- Introducer Cannula
- 1- MIC-KEY* ENFit™ Extension Set with Right Angle Connector and double male ENFit™ Connector, 30cm length
- 1- MIC-KEY* ENFit™ Bolus Extension Set with Straight Connector and male ENFit™ connector, 30cm length
- 1- 6ml Luer-Slip™ Syringe
- 1- 35ml female ENFit™ syringe
- 4- Gauze Pads

HALYARD* MIC-KEY* J-TUBE


The **MIC-KEY* J-TUBE** is indicated for patients requiring jejunal feeding only, but who want an inconspicuous external segment. With an external bolster designed like the widely prescribed MIC-KEY* G-Tube.

- Medical-Grade Silicone Construction
- Low-Profile Design
- Inflatable Silicone Internal Retention Balloon
- High performing Anti-Reflux Valve
- SECUR-LOK* Extension Set Connector Mechanism with ENFit™ Connectors
- Trimmable Distal Tip
- Radiopaque Striping
- Not formulated with Natural Rubber Latex nor DEHP
- Gamma Sterilised
- 1 Each/Packaging

HALYARD* MIC-KEY* J-TUBES

| Size (Outer Dia.) (X) | Stoma Length (Y) | Product Code |
|-----------------------|------------------|--------------|
| 14 Fr | From 0.8 to 4.5 | 8230-X-Y |
| 18 Fr | From 0.8 to 4.5 | 8230-X-Y |

Available stoma lengths (cm): 0.8 /1.0 /1.2 /1.5 /1.7 /2.0 /2.3 /2.5 /2.7 /3.0 /3.5 /4.0 /4.5




HALYARD* MIC* J-TUBE

The **MIC* J-TUBE** is indicated for patients requiring jejunal feeding only. The distal tip can be trimmed to suit individual patient needs.

- Medical-Grade Silicone Construction
- Inflatable Silicone Internal Retention Balloon
- SECUR-LOK* External Retention Ring
- Male ENFit™ Feeding Port
- Trimmable Distal Tip
- Radiopaque Striping
- Not formulated with Natural Rubber Latex nor DEHP
- Gamma Sterilised
- 2 Eaches/Packaging

HALYARD* MIC* J-TUBES

| Size (Outer Dia.) | Balloon Volume | Product Code |
|-------------------|----------------|--------------|
| 12 Fr | 3-5ml Balloon | 8200-12LV |
| 14 Fr | 7-10ml Balloon | 8200-14 |
| 16 Fr | 7-10ml Balloon | 8200-16 |
| 18 Fr | 7-10ml Balloon | 8200-18 |
| 20 Fr | 7-10ml Balloon | 8200-20 |
| 22 Fr | 7-10ml Balloon | 8200-22 |
| 24 Fr | 7-10ml Balloon | 8200-24 |



HALYARD* MIC* JEJUNOSTOMY TUBE

The **MIC* JEJUNOSTOMY TUBE** is surgically placed directly into the small bowel, utilising the Witzel Tunnel technique for secure placement and to help minimise leakage. Featuring a large French 14 tubing to help prevent clogging.

- Medical-Grade Silicone Construction
- Suture Wings
- French 14 Tubing (Large 9 French Internal Diameter)
- Radiopaque Stripe
- Gamma Sterilised
- Male Luer Lock Feeding Connector
- 1 Each/Packaging

Kit includes:

- 1- Jejunostomy Tube
- 1- Subcutaneous Extraducer
- 1- Feeding Port Adapter
- 1- Luer-Lok® Cap

HALYARD* MIC* JEJUNOSTOMY TUBES

| Size (Outer Dia.) | Product Code |
|-------------------|--------------|
| 14 Fr | 0301-14 |

GASTRO-JEJUNAL ACCESS FEEDING TUBES

HALYARD* Gastro-Jejunal tubes allow simultaneous enteral nutrition in the small intestine and gastric decompression. These tubes are available in low-profile and conventional, standard-length designs.



HALYARD* MIC-KEY* GJ-TUBE

The unique **MIC-KEY* GJ-TUBE** is designed for patients who require simultaneous gastric decompression and jejunal feeding. Available in several sizes and jejunal lengths to meet individual patient needs.

- Medical-Grade Silicone Construction
- Low-Profile Design
- Tapered Distal Tip matching tube size
- Inflatable Silicone Internal Retention Balloon (3-5ml, tubes with 45 cm Jejunal length have a 7-10ml balloon)
- Dual Ports: Male ENFit™ Gastric Decompression and Jejunal Feeding Port
- Multiple Gastric & Jejunal Exit Ports
- Tungsten Weighted Jejunal Portion
- Radiopaque Striping
- Gamma sterilised
- 1 Each/Packaging

HALYARD* MIC-KEY* GJ-TUBES

| Size (Outer Dia.) | Stoma Length | Jejunal Length 15 CM | 22 CM | 30 CM | 45 CM |
|-------------------|--------------|----------------------|----------------|----------------|----------------|
| 16 Fr | 1.0cm | 8270-16-1.0-15 | 8270-16-1.0-22 | | |
| 16 Fr | 1.2cm | 8270-16-1.2-15 | 8270-16-1.2-22 | 8270-16-1.2-30 | 8270-16-1.5-45 |
| 16 Fr | 1.5cm | 8270-16-1.5-15 | 8270-16-1.5-22 | 8270-16-1.5-30 | 8270-16-1.5-45 |
| 16 Fr | 1.7cm | 8270-16-1.7-15 | 8270-16-2.0-22 | 8270-16-1.7-30 | 8270-16-1.7-45 |
| 16 Fr | 2.0cm | | | 8270-16-2.0-30 | 8270-16-2.0-45 |
| 16 Fr | 2.3cm | | | 8270-16-2.3-30 | 8270-16-2.3-45 |
| 16 Fr | 2.5cm | | | 8270-16-2.5-30 | 8270-16-2.5-45 |
| 16 Fr | 2.7cm | | | 8270-16-2.7-30 | 8270-16-2.7-45 |
| 16 Fr | 3.0cm | | | | 8270-16-3.0-45 |
| 16 Fr | 3.5cm | | | | 8270-16-3.5-45 |
| 18 Fr | 1.2cm | | 8270-18-1.2-22 | 8270-18-1.2-30 | |
| 18 Fr | 1.5cm | | 8270-18-1.5-22 | 8270-16-1.5-30 | 8270-18-1.5-45 |
| 18 Fr | 1.7cm | | 8270-18-1.7-22 | 8270-18-1.7-30 | 8270-18-1.7-45 |
| 18 Fr | 2.0cm | | 8270-18-2.0-22 | 8270-18-2.0-30 | 8270-18-2.0-45 |
| 18 Fr | 2.3cm | | | 8270-18-2.3-30 | 8270-18-2.3-45 |
| 18 Fr | 2.5cm | | | 8270-18-2.5-30 | 8270-18-2.5-45 |
| 18 Fr | 2.7cm | | | 8270-18-2.7-30 | 8270-18-2.7-45 |
| 18 Fr | 3.0cm | | | | 8270-18-3.0-45 |
| 18 Fr | 3.5cm | | | | 8270-18-3.5-45 |
| 22 Fr | 1.5cm | | | | 8270-22-1.5-45 |
| 22 Fr | 1.7cm | | | | 8270-22-1.7-45 |
| 22 Fr | 2.0cm | | | | 8270-22-2.0-45 |
| 22 Fr | 2.3cm | | | | 8270-22-2.3-45 |
| 22 Fr | 2.5cm | | | | 8270-22-2.5-45 |
| 22 Fr | 2.7cm | | | | 8270-22-2.7-45 |
| 22 Fr | 3.0cm | | | | 8270-22-3.0-45 |
| 22 Fr | 3.5cm | | | | 8270-22-3.5-45 |

Kit Includes:

- 1- MIC-KEY* Low-Profile GJ-tube
- 1- Introducer Cannula
- 1- MIC-KEY* ENFit™ Extension Set with Right Angle Connector and double male ENFit™ Connector, 30cm length
- 1- MIC-KEY* ENFit™ Bolus Extension Set with Straight Connector and male ENFit™ connector, 30cm length
- 1- 6ml Luer-Slip™ Syringe
- 1- 35ml female ENFit™ syringe
- 4- Gauze Pads




HALYARD* MIC* GJ-TUBE (Endoscopic/Fluoroscopic Placement)

The unique **MIC* GASTRO-JEJUNAL TUBE** is also designed for patients who require simultaneous gastric decompression and jejunal feeding. Available in several sizes and jejunal lengths to meet individual patient needs.

- Medical-Grade Silicone Construction
- Inflatable Silicone Internal Retention Balloon (7-10ml)
- SECUR-LOK* External Retention Ring
- Dual Ports: Male ENFit™ Gastric Decompression and Jejunal Feeding Port
- Multiple Gastric & Jejunal Exit Ports
- Tungsten Weighted Jejunal Portion
- Tapered Distal Tip matching tube size
- Radiopaque Striping
- Gamma Sterilised
- 1 Each/Packaging

HALYARD* MIC* GJ-TUBES (Endoscopic/Fluoroscopic Placement)

| Size (Outer Dia.) | Jejunal Length 15 CM | 22 CM | 30 CM | 45 CM |
|-------------------|----------------------|------------|------------|---------|
| 16 Fr | 8650-16-15 | 8650-16-22 | 8650-16-30 | 8650-16 |
| 18 Fr | | 8650-18-22 | 8650-18-30 | 8650-18 |
| 22 Fr | | | | 8650-22 |



HALYARD* MIC* GJ-TUBE (Surgical Placement)


The unique **MIC* GJ-TUBE** is also designed for patients who require simultaneous gastric decompression and jejunal feeding. Tube for surgical placement, available in several French sizes and 45 cm jejunal length.

- Medical-Grade Silicone Construction
- Inflatable Silicone Internal Retention Balloon (7-10 ml)
- SECUR-LOK* External Retention Ring
- Dual Ports: Male ENFit™ Gastric Decompression and Jejunal Feeding Port
- Multiple Gastric & Jejunal; Exit Ports
- Tungsten Weighted Jejunal Portion
- Tapered Distal Tip matching tube size
- Radiopaque Striping
- ETO Sterilised
- 1 Each/Packaging

Kit includes:
 1- MIC* GJ-tube
 1- Abdominal Trocar and Sheath
 1- Slotted Transpyloric Cannula
 1- 5x5cm Slit Gauze Sponge

HALYARD* MIC* GJ-TUBES (Surgical Placement)

| Size (Outer Dia.) | Jejunal Length 45 CM |
|-------------------|----------------------|
| 16 Fr | 8660-16 |
| 18 Fr | 8660-18 |
| 22 Fr | 8660-22 |



HALYARD* MIC* GASTRO-ENTERIC TUBE

The **MIC* GASTRO-ENTERIC TUBE** is a precursor to the very popular MIC* GJ-Tube. Also designed for simultaneous gastric decompression and jejunal feeding and for patients designated to a smaller distal jejunal segment.

- Medical-Grade Silicone Construction
- Inflation Silicone Internal Retention Balloon (7-10ml, tubes with 25.4cm Jejunal length have a 3-5 ml balloon)
- SECUR-LOK* External Retention Ring
- Tungsten weighted jejunal tube
- Dual ports: Male ENFit™ Gastric Decompression and Jejunal Feeding Port
- Multiple Gastric & Jejunal Exit Ports
- Radiopaque Striping and full opaque jejunal tubing
- Gamma Sterilised
- 2 Eaches/Packaging

HALYARD* MIC* GASTRO-ENTERIC TUBES

| Size (Outer Dia.) | 25.4 CM | 57.9 CM |
|-------------------|-----------|---------|
| 16 Fr | 8210-16LV | 8210-16 |
| 18 Fr | 8210-18LV | 8210-18 |
| 20 Fr | | 8210-20 |
| 22 Fr | | 8210-22 |
| 24 Fr | | 8210-24 |
| 26 Fr | | 8210-26 |
| 28 Fr | | 8210-28 |
| 30 Fr | | 8210-30 |



HALYARD* INTRODUCER KIT FOR J- AND GJ-TUBE PLACEMENT

Kits include tools and devices to facilitate primary endoscopic, fluroscopic or laparoscopic placement of balloon retained Jejunostomy and Gastro-Jejunostomy tubes such as MIC* J-/GJ-Tubes and MIC-KEY* J-/GJ-tubes.

- ETO Sterilised
- 1 Each/Packaging

HALYARD* INTRODUCER KITS FOR J- AND GJ-TUBE PLACEMENT

| Description | Product Code |
|--|--------------|
| Introducer Kit for J- and GJ-Tube placement, with 20 Fr Dilator (for use with MIC* J-/GJ-Tube Fr 16 and MIC-KEY* J-/GJ-Tube Fr 16) | 98437 |
| Introducer Kit for J- and GJ-Tube placement, with 22 Fr Dilator (for use with MIC* J-/GJ-Tube Fr 18 and MIC-KEY* J-/GJ-Tube Fr 18) | 98438 |

Kit includes:

- 1- Gastrointestinal Anchor Set with 4 Pre-loaded SAFE-T-PEXY* T-Fasteners
- 1- Haemostat
- 1- Safety Introducer Needle
- 1- Seeking Catheter Fr 5
- 1- Safety Scalpel #11
- 1- Serial Dilator with Peel-away Sheath
- 2- Syringes 12ml
- 1- Stoma Measuring Device

Feeding tube not included in kit.

ACCESSORIES

MIC*/MIC-KEY* Accessories support gastric, jejunal and gastro-jejunal access feeding tubes.



- SECUR-LOK* "Key" Connection Mechanism
- Double male ENFit™ Connector (except MIC-KEY* Bolus Extension Set)
- Individually Packaged Extension Sets
- Tubing Clamp (except on medication set)
- Not formulated with Natural Rubber Latex nor DEHP
- Non Sterile

HALYARD* MIC-KEY* EXTENSION SETS

For use with the MIC-KEY* G / J / GJ-Tubes

| Description | Length | Packaging | Product Code |
|---|--------|-----------|--------------|
| MIC-KEY* Extension Set ENFit™ with Right-Angle Connector | 30cm | 5 Each | 0141-12 |
| MIC-KEY* Extension Set ENFit™ with Right-Angle Connector | 60cm | 5 Each | 0141-24 |
| MIC-KEY* Medication Set ENFit™ with Right-Angle Connector | 5cm | 5 Each | 0142-02 |
| MIC-KEY* Bolus Extension Set ENFit™ with Straight Connector | 30cm | 5 Each | 0143-12 |
| MIC-KEY* Bolus Extension Set ENFit™ with Straight Connector | 60cm | 5 Each | 0143-24 |
| MIC-KEY* Stoma Measuring Device | | 10 Each | 98460 |

HALYARD* MIC* JEJUNOSTOMY FEEDING TUBE ACCESSORIES

For use with the MIC* Jejunostomy Feeding Tube

| Description | Packaging | Product Code |
|---------------------------------------|-----------|--------------|
| MIC* Jejunostomy Tube Repair Kit | 1 Each | 0301-50 |
| MIC* Jejunostomy Female Luer-Lok Caps | 5 Each | 0301-55 |
| MIC* Jejunostomy Feeding Port Adapter | 5 Each | 0301-60 |

HALYARD* MIC* PEG FEEDING CONNECTORS

For use with MIC* Percutaneous Endoscopic Gastrostomy Feeding Tubes

- Not formulated with Natural Rubber Latex nor DEHP
- Non Sterile
- 1 Each/Packaging

| Description | Product Code |
|-------------------------------------|--------------|
| Male ENFit™ Feeding Connector Fr 14 | 8135-14 |
| Male ENFit™ Feeding Connector Fr 20 | 8135-20 |
| Male ENFit™ Feeding Connector Fr 24 | 8135-24 |



KNOWLEDGE NETWORK* Clinical Education
Knowledgeable Customer Support
Expert Sales Force
Tools & Best Practices
Clinical Research
Commitment to Excellence

For more information, please send an email to
customerservice.uk.ie@hyh.com or visit
www.halyardhealth.co.uk



ENFit™ : Registered Trademark of Global Enteral Device Supplier Association, Inc.
Used with their permission.

*Registered Trademark or Trademark of Halyard Health, Inc. or its affiliates.
© 2017 HYH. All rights reserved. HC113-00-UK

www.halyardhealth.co.uk