



## Polyisoprene Surgical Gloves, Powder Free

Our GloveOn® Victor polyisoprene surgical glove provides surgeons with a fine touch for surgical procedures requiring a high tactile response. Combined with great dexterity, comfort and feel, GloveOn® Victor allows you to concentrate fully on the job at hand. Best used when worn in conjunction with GloveOn® Victor Under Glove for a complete double donning solution that aids in puncture detection.



Tested in accordance with EN ISO 374-5:2016  
Resistance to Bacteria and Fungi - pass  
Resistance to Viruses - pass



Protection against particulate radioactive contamination (Excluding Clause 4.3)

Physical Dimensions				
Length (mm)	290 ± 10			
Cuff Thickness (mm)	0.15 ± 0.04			
Palm Thickness (mm)	0.19 ± 0.04			
Finger Thickness (mm)	0.21 ± 0.04			
Physical Properties	Before Aging	After Aging		
	AS/NZS 4179	Typical Value	AS/NZS 4179	Typical Value
Force at break (N)	≥ 12.5	16	≥ 9.5	14
Elongation at break (%)	≥ 700	1053	≥ 550	1069
Modulus at 300% elongation (N)	≤ 2.0	1.2	-	-
Inspection Levels & AQL	Inspection Level	AQL		
Watertightness	G1	0.65		
Physical dimensions	S2	4.0		
Physical properties	S2	4.0		
Visual inspection (major)	S4	2.5		
Visual inspection (minor)	S4	4.0		
Particulate residue	N = 5	≤ 2 mg/glove		
Seal strength	N = 20	≥ 1.2 N		
Dye penetration test	N = 20	Accept 0 Reject 1		

### REORDER CODE

SPI2855	Size 5.5
SPI2860	Size 6.0
SPI2865	Size 6.5
SPI2870	Size 7.0
SPI2875	Size 7.5
SPI2880	Size 8.0
SPI2885	Size 8.5
SPI2890	Size 9.0

### FEATURES

- Powder free • Not made with natural rubber latex • Textured • Curved fingers
- Hand specific • Natural white colour
- Sterilised using gamma irradiation

### PACKAGING

1 pair per pouch  
50 pairs per box  
4 boxes per carton

### REGULATORY COMPLIANCE

ARTG 291855, FDA 510 (k), MDD 93/42/EEC

### STANDARDS

AS/NZS 4179, ASTM D3577, ASTM D5151, ASTM D6124, ASTM D6978, ASTM D7102, ASTM F1671, EN 455 part 1, 2 & 3, EN 16523-1, ISO 10282, ISO 10993 part 5 & 10

### MANUFACTURING ACCREDITATIONS

ISO 9001, ISO 13485, EN ISO 13485

Measured breakthrough time (minutes)	>10	>30	>60	>120	>240	>480
Permeation performance level	1	2	3	4	5	6

Chemical	EN 16523-1:2015 Permeation Level	EN 374-4:2013 Mean Degradation (%)
I Ethyl Acetate	6	49.4
K 40% Sodium Hydroxide	6	-0.9
M 65% Nitric Acid	5	25.1
P 30% Hydrogen Peroxide	6	-8.0
T 37% Formaldehyde	6	-31.0

A brand by



**Mun Australia**  
Toll free: 1800 456 837  
munglobal.com.au



Proud Supporter of  
Birthing Kit Foundation (Australia)



Protection Always On