



SAFETOUCH™ NEEDLE

STERILE, SINGLE-USE WITH SAFETY DEVICE
FOR INJECTION OR ASPIRATION



SAFETOUCH™ NEEDLE

STERILE, SINGLE-USE WITH SAFETY DEVICE FOR INJECTION OR ASPIRATION

- Transparent safety mechanism activated by one hand
- Audible click confirms locking of safety device
- Does not require a change in puncture technique
- Ultra sharp three-beveled needle to minimize discomfort
- Siliconization of needle minimizes penetration and gliding force
- Color-coded to conform to ISO 6009
- Not made with latex, DEHP or PVC
- Available in 16-27 G
- Available in thin or regular wall
- Available with needle cap, short (42 mm) or long (54 mm)
- For use by healthcare professionals only



Class IIa Medical Devices
Rule 6
MDD 93/42/EEC
UMDNS 12729

PRODUCT COMPLIANCE

- CE-marked, Class IIa Medical Devices, Rule 6, MDD 93/42/EEC, UMDNS 12729
- Complies with the following norms, directives, and regulations:
 - ISO 80369-7:2016
 - ISO 7864:1993
 - ISO 11138-1:2006
 - ISO 2859-1:1999
 - ISO 7886-1:1993
 - ISO 23908:2011
 - ISO 6009:1992
 - ISO 10993-10:2010
- Labels contain 26 languages

MANUFACTURING DETAILS

Legal manufacturer: Nipro Corporation

Country of origin: Thailand

STERILIZATION AND SHELF LIFE

Sterilization

EtO (Ethylene oxide)



Single-use only



Indicator sticker


Each outer box (or shipping carton) contains a chemical indicator that indicates sterility. This indicator comes in the form of a [violet] sticker that turns [blue] when sterilized.

Shelf life

5 years



MATERIALS USED



**Not made with
Latex, DEHP or PVC**

PRODUCT RANGE OVERVIEW

Reference code	Description	Gauge (G)	Length (mm)	Color
HN-SF-1625-ET	SAFETOUCH Needle 16 G x 1" (1.2x25 mm) LR TW	16	25	Grey
HN-SF-1825-ET	SAFETOUCH Needle 18 G x 1" (1.2x25 mm) LR TW	18	25	Pink
HN-SF-1832-ET	SAFETOUCH Needle 18 G x 1 1/4" (1.2x32 mm) LR TW	18	32	Pink
HN-SF-1838-ET	SAFETOUCH Needle 18 G x 1 1/2" (1.2x38 mm) LR TW	18	38	Pink
HN-SF-1925-ET	SAFETOUCH Needle 19 G x 1" (1.1x25 mm) LR TW	19	25	Brown
HN-SF-1932-ET	SAFETOUCH Needle 19 G x 1 1/4" (1.1x32 mm) LR TW	19	32	Brown
HN-SF-1938-ET	SAFETOUCH Needle 19 G x 1 1/2" (1.1x38 mm) LR TW	19	38	Brown
HN-SF-2025-ET	SAFETOUCH Needle 20 G x 1" (0.9x25 mm) LR TW	20	25	Yellow
HN-SF-2038-ET	SAFETOUCH Needle 20 G x 1 1/2" (0.9x38 mm) LR TW	20	38	Yellow
HN-SF-2116-ET	SAFETOUCH Needle 21 G x 5/8" (0.8x16 mm) LR TW	21	16	Green
HN-SF-2125-ET	SAFETOUCH Needle 21 G x 1" (0.8x25 mm) LR TW	21	25	Green
HN-SF-2132-ET	SAFETOUCH Needle 21 G x 1 1/4" (0.8x32 mm) LR TW	21	32	Green
HN-SF-2138-ET	SAFETOUCH Needle 21 G x 1 1/2" (0.8x38 mm) LR TW	21	38	Green
HN-SF-2225-ET	SAFETOUCH Needle 22 G x 1" (0.7x25 mm) LR TW	22	25	Black
HN-SF-2232-ET	SAFETOUCH Needle 22 G x 1 1/4" (0.7x32 mm) LR TW	22	32	Black
HN-SF-2238-ET	SAFETOUCH Needle 22 G x 1 1/2" (0.7x38 mm) LR TW	22	38	Black
HN-SF-2316-ET	SAFETOUCH Needle 23 G x 5/8" (0.6x16 mm) LR TW	23	16	Blue
HN-SF-2325-ET	SAFETOUCH Needle 23 G x 1" (0.6x25 mm) LR TW	23	25	Blue
HN-SF-2332-ET	SAFETOUCH Needle 23 G x 1 1/4" (0.6x32 mm) LR TW	23	32	Blue
HN-SF-2338-ET	SAFETOUCH Needle 23 G x 1 1/2" (0.6x38 mm) LR RW	23	38	Blue
HN-SF-2425-ET	SAFETOUCH Needle 24 G x 1" (0.55x25 mm) LR TW	24	25	Purple
HN-SF-2432-ET	SAFETOUCH Needle 24 G x 1 1/4" (0.55x32 mm) LR TW	24	32	Purple
HN-SF-2513-ET	SAFETOUCH Needle 25 G x 1/2" (0.5x13 mm) LR TW	25	13	Orange
HN-SF-2516-ET	SAFETOUCH Needle 25 G x 5/8" (0.5x16 mm) LR TW	25	16	Orange
HN-SF-2525-ET	SAFETOUCH Needle 25 G x 1" (0.5x25 mm) LR TW	25	25	Orange

Reference code	Description	Gauge (G)	Length (mm)	Color
HN-SF-2532-ET	SAFETOUCH Needle 25 G x 1 1/4" (0.5x32 mm) LR TW	25	32	■
HN-SF-2613-ET	SAFETOUCH Needle 26 G x 1/2" (0.45x13 mm) LR RW	26	13	■
HN-SF-2622-ET	SAFETOUCH Needle 26 G x 7/8" (0.45x22 mm) LR RW	26	22	■
HN-SF-2625-ET	SAFETOUCH Needle 26 G x 1" (0.45x25 mm) LR RW	26	25	■
HN-SF-2713-ET	SAFETOUCH Needle 27 G x 1/2" (0.4x13 mm) LR RW	27	13	■
HN-SF-2719-ET	SAFETOUCH Needle 27 G x 3/4" (0.4x19 mm) LR RW	27	19	■
HN-SF-2732-ET	SAFETOUCH Needle 27 G x 1 1/4" (0.4x32 mm) LR RW	27	32	■
HN-SF-2738-ET	SAFETOUCH Needle 27 G x 1 1/2" (0.4x38 mm) LR RW	27	38	■

PACKAGING DETAILS

Reference code	Inner box				Outer box				
	Quantity	L (mm)	W (mm)	H (mm)	Quantity	L (mm)	W (mm)	H (mm)	Weight (kg)
HN-SF-XXXX-ET	50	255	112	63	500	342	246	287	2,800

Languages

26 languages on IFU, inner box, and outer box:

English (EN), Bulgarian (BG), Czech (CS), Danish (DA), German (DE), Greek (EL), Spanish (ES), Estonian (ET), Finnish (FI), French (FR), Croatian (HR), Hungarian (HU), Italian (IT), Lithuanian (LT), Latvian (LV), Dutch (NL), Norwegian (NO), Polish (PL), Portuguese (PT), Romanian (RO), Russian (RU), Slovakian (SK), Slovenian (SL), Serbian (SR), Swedish (SV), Turkish (TR)

Transport conditions

Closed and dry

Storage conditions

Open the packaging only immediately before use to guarantee sterility. Store in a dry place away from direct heat and sunlight.

Label details

	Individual unit	Inner box	Outer box
Product name	✓	✓	✓
Product size	✓	✓	✓
Quantity (pieces)	—	✓	✓
REF	✓	✓	✓
LOT	✓	✓	✓
	✓	✓	✓
STERILE EO	✓	✓	✓
CE	✓	✓	✓
	✓	✓	✓

	Individual unit	Inner box	Outer box
	✓	✓	✓
	✓	✓	✓
	✓	✓	✓
EC REP	✓	✓	✓
	—	✓	✓
	✓	✓	✓
	✓	✓	✓
	✓	✓	✓

