

SureSafe Surgical Headwear



Bouffant Caps



Theatre Cap

NEW



Headscarf

NEW

Features

Single layer spunbond polypropylene fabric

Provides minimum fluid resistance

Soft and comfortable

Latex free, lint-free and anti-static

Non-sterile

PRODUCT CODE	PRODUCT DESCRIPTION	SIZE	COLOUR	SKU/TU (MOQ)	SKU/CTN
PA31021B	SureSafe Bouffant Cap	21" (53CM)	Blue	500	500
PA31021R	SureSafe Bouffant Cap (Use SM 10125114)	21" (53CM)	Red	500	500
PA31021W	SureSafe Bouffant Cap	21" (53CM)	White	500	500
PA37002B	SureSafe Headscarf	96 x 96cm	Blue	250	250
PA37002P	SureSafe Headscarf	96 x 96cm	Pink	250	250
PA37001B	SureSafe Theatre Cap, Tie Back	OSFA	Blue	250	250
PA37001P	SureSafe Theatre Cap, Tie back	OSFA	Pink	250	250