



# ENDOSCOPY PRODUCTS AND CONSUMABLES



ESSENTIAL TO HEALTH

# Drawer-Style Trap



## Product Description

Smooth, intuitive, and user-friendly, the Drawer-Style Trap facilitates the processing of unlimited specimens and offers superior visualisation to ensure secure recovery.

EBOS Code	Supplier Code	Description	UOM
10022300	VIT-SUEPTP-02D	ENDO-SAFIER DISPOSABLE POLYP TRAP WITH 2 DRAWERS	Box/10

# 5-Chamber Trap



## Product Description

The 5-Chamber Trap has 5 individually removable and numbered chambers. Each trap uses a double layer of fine mesh, ensuring polyps are not lost and resulting in superior specimen retention and separation from waste fluids.

EBOS Code	Supplier Code	Description	UOM
10022301	VIT-SUEPTP-05R	DISPOSABLE POLYP TRAP, 5 REMOVEABLE CHAMBERS	Box/10

# Single Trap



## Product Description

The Single Trap confidently catches your specimen when only one is needed. Simple to use and great value for money.

EBOS Code	Supplier Code	Description	UOM
10022299	VIT-SUEPTP-001	ENDO-SAFIER DISPOSABLE POLYP TRAP SINGLE	Box/25

# Single Use Proctoscope with Light Source

## Key Features

- Extra-large diameter (30mm) allowing more room for surgical instruments.
- Fenestrated to allow easier access to mucosa.
- LED light source.
- 32 mm lumen to accommodate most devices.




EBOS Code	Supplier Code	Description	UOM
10022297	VIT-SUE-PRO-LIT	PROCTOSCOPE SNG USE WTH LIGHT SOURCE	Box/25

# Precision Clip Endoscopic Hemoclip

## Product Description

The new Precision Hemoclip combines robust holding strength with a new level of clip placement control.

This clip offers genuine one-to-one rotation, controllable at the port and performs in tight loops. The clip can be repositioned repeatedly prior to deployment, which is simple and instinctive without requiring excessive force.



## Key Features

- One-to-one rotation, even in tight loops.
- Repositionable - clip can be opened and closed multiple times prior to final deployment.
- Handle locks after deployment leaving the tip smooth (without a hook), avoiding damage to scope or tissue after releasing clips.
- Release system is sensitive and intuitive and does not require excessive force to effect clip deployment.

EBOS Code	Supplier Code	Description	UOM
10022293	VIT-SUE-HCL-016	HEMOCLIP 16MM JAW	Box/5
10022292	VIT-SUE-HCL-011	HEMOCLIP 11MM JAW	Box/5

## Specifications

- Working Length 2300mm
- Required Channel  $\geq$  2.8mm diameter
- Jaw Opening 11mm or 16mm
- Units per Box 5 units
- Rebate Code VM002



# Precision Snare


## Key Features

- 360°1:1 rotation makes operation more flexible even in tight loops.
- Precision sharpness for a clean cut, resulting in no ragged edges or tears.
- Ergonomic handle for smooth operation requiring minimal effort.

## Specifications

- Working Length 2300mm
- Required Channel  $\geq$  2.8mm diameter
- Loop Size 10 mm
- Units per Box 10

EBOS Code	Supplier Code	Description	UOM
10022298	VIT-SUE-SPO-010	PRECISION SNARE, 10MM LOOP, DIAMOND SHAPE, 2300MM LENGTH	Box/10



# Snares

EBOS Code	Supplier Code	Description	UOM
10022303	VIT-SUE-PSO-015	SNARE 15MM	Box/10
10022304	VIT-SUE-PSO-030	SNARE 30MM	Box/10

# Vitramed Cleaning Brushes

## Product Description

Vitramed cleaning brushes offer reliability and consistent performance and are robust enough to resist bending and/or buckling under pressure yet still represent great value. They can be packed individually or made into kits upon request.

### Long Brush (Double Ended)

- Catheter Length 2250mm
- Brush Diameter 6mm
- Brush Diameter 6mm

### Valve Brush (Double Ended)

- Overall Length 150mm
- Brush Diameter (Large End) 10mm
- Brush Diameter (Small End) 5mm



EBOS Code	Supplier Code	Description	UOM
10022266	VIT-SUE-PSO-015	CLEANING BRUSH KIT	Box/50
10022267	VIT-SUE-PSO-030	CLEANING BRUSH KIT (INC TOOTHBRUSH)	Box/40

# Vitramed Bougie Dilators

## Product Description

Vitramed offers single use oesophageal dilators in a range of sizes from 5mm through to 20mm in diameter. Our bougies are clearly marked at 1cm intervals for optimum accuracy and carry a radio-opaque marker. Vitramed bougies are a low-cost alternative to reusable oesophageal dilators, removing the hassle and expense of reprocessing.



EBOS Code	Supplier Code	Description	UOM
10022274	VIT-SUE-DDL-005	BOUGIE DILATOR 5MM	Each
10022275	VIT-SUE-DDL-006	BOUGIE DILATOR 6MM	Each
10022276	VIT-SUE-DDL-007	BOUGIE DILATOR 7MM	Each
10022277	VIT-SUE-DDL-008	BOUGIE DILATOR 8MM	Each
10022278	VIT-SUE-DDL-009	BOUGIE DILATOR 9MM	Each
10022279	VIT-SUE-DDL-010	BOUGIE DILATOR 10MM	Each
10022280	VIT-SUE-DDL-011	BOUGIE DILATOR 11MM	Each
10022281	VIT-SUE-DDL-012	BOUGIE DILATOR 12MM	Each

EBOS Code	Supplier Code	Description	UOM
10022282	VIT-SUE-DDL-013	BOUGIE DILATOR 13MM	Each
10022283	VIT-SUE-DDL-014	BOUGIE DILATOR 14MM	Each
10022284	VIT-SUE-DDL-015	BOUGIE DILATOR 15MM	Each
10022285	VIT-SUE-DDL-016	BOUGIE DILATOR 16MM	Each
10022286	VIT-SUE-DDL-017	BOUGIE DILATOR 17MM	Each
10022287	VIT-SUE-DDL-018	BOUGIE DILATOR 18MM	Each
10022288	VIT-SUE-DDL-019	BOUGIE DILATOR 19MM	Each
10022289	VIT-SUE-DDL-020	BOUGIE DILATOR 20MM	Each

# Single Use Guide Wire for Bougie Dilators

## Product Description

Guide wires are an accompaniment to the bougies. Made of stainless steel and carry markings every 10cm.

EBOS Code	Supplier Code	Description	UOM
10022290	VIT-SUE-GDW-000	GUIDEWIRE 0.9MM DIAMETER, 2000MM LENGTH	Box/5
10022291	VIT-SUE-GDW-002	GUIDEWIRE 0.9MM DIAMETER, 2400MM LENGTH	Box/5



# Cold/Hot Biopsy Forceps

## Product Description


Vitramed offer a range of single-use, cold and hot biopsy forceps to suit all budgets and preferences. 4.8mm<sup>3</sup> fenestrated smooth cup basket with compact quad pole jaw structure. These forceps are laser cut to ensure a clean bite every time.

### Key Features (Cold)

- Larger 7.8mm<sup>3</sup> alligator cup.
- Standard catheter length is 2300mm but other lengths and diameters are available by request.

### Key Features (Hot)

- DuPont Teflon PTFE cable coating for easiest insertion through the channel.
- Extreme resistance to abrasions, acid or alkaline environments.



EBOS Code	Supplier Code	Description	UOM
10022268	VIT-SUE-CFP-001	COLD FORCEP STANDARD SMOOTH CUP 2.3MM DIAMETER, 2300MM LENGTH	Box/10
10022271	VIT-SUE-CFP-LG1	COLD FORCEP ALLIGATOR JAW 2.3MM DIAMETER, 2300MM LENGTH	Box/10
10022273	VIT-SUE-CFP-M01	COLD FORCEP (PAEDIATRIC) 1.8MM DIAMETER, 1600MM LENGTH	Box/10
10022269	VIT-SUE-GFC-230	GRASPING FORCEP	Box/10
10022270	VIT-SUE-CFP-002	COLD FORCEP SMOOTH CUP AGS	Box/10
10022272	VIT-SUE-CFP-LG2	COLD FORCEP ALLIGATOR JAW AGS	Box/10
10022302	VIT-SUE-HFP-000	HOT BIOPSY FORCEPS 230CM X 2.3MM	Box/10

# Button and Valve Set

## Product Description

Vitramed single-use button and valve set is presented in a box of 20 and represent the right combination of price and performance. Tested to ensure a tight seal and to minimise leakage. The set contains water, suction, biopsy, water jet adapter and cleaning valve.

EBOS Code	Supplier Code	Description	UOM
10022296	VIT-SUE-OCB-005	OLYMPUS BUTTON & VALVE SET	Box/20 Sets



# Accessories

EBOS Code	Supplier Code	Description	UOM
10022264	VIT-SUE-CFP-001	HELIPROBE BREATH CARDS	Box/5
10022265	VIT-SUE-CFP-LG1	HELICAP	Box/10
10022306	VIT-SUE-CFP-M01	RETRIVAL NET 30MM	Box/10
10022294	VIT-SUE-GFC-230	VITRAMED INJ NEEDLE 5MM 23G	Box/10
10022295	VIT-SUE-CFP-002	VITRAMED INJ NEEDLE 5MM 25G	Box/10





EBOS HEALTHCARE:

300 Victoria Street,  
Wetherill Park,  
NSW 2164

Free Phone: 1800 269 534

Free Fax: 1800 810 257

Email: [customerservice@ebosgroup.com.au](mailto:customerservice@ebosgroup.com.au)

Online Ordering: [www.eboshealthcare.com.au](http://www.eboshealthcare.com.au)

WAREHOUSE LOCATIONS:

New South Wales:

300 Victoria Street,  
Wetherill Park,  
NSW 2164

Queensland:

Unit 2, 39 Distribution Street  
Motorway Business Park  
Larapinta  
QLD 4110

Victoria:

34-78 South Park Drive,  
Dandenong South,  
VIC 3175

Western Australia:

2/85 Sheffield Road  
Welshpool  
WA 6106

South Australia:

75-83 Hardys Road  
Underdale  
SA 5032

Tasmania:

7A Pear Avenue  
Derwent Park  
TAS 7009



ESSENTIAL TO HEALTH