



Water Resistant



Contours to Body



Transparent Film

A transparent adhesive polyurethane film with moisture vapour permeability which adheres to the skin and forms a flexible and comfortable layer. AsGuard Clear conforms to the skin and allows easy visibility to wound healing.

CUT TO LENGTH

Measured grid backing allows cut to correct size as required

TRANSPARENT FILM

Allows visualisation of covered tissue

SHOWER AND INFECTION RESISTANT

Impermeable to help protect site against infection and allows showering with the dressing in-situ



Suitable for:

Wounds with little or no exudate

Protection and securement of IV lines

Dressing reinforcement

Reorder	Size	Pcs/Unit	Units/Ctn
ISB010	5cm x 10m	1	24
ISB011	10cm x 10m	1	12

